

## SYNCERITY™ 1024 × 256

### Open-Electrode TE-Cooled CCD Detector for UV-VIS-NIR

#### Lowest Noise & Highest Dynamic Range in its Class

#### Research Grade OEM Spectroscopy Camera for 200 nm–1050 nm

Flexible pricing structure scalable with performance and OEM volume requirements

The deep cooled Open Electrode 1024 × 256 CCD Camera combines affordability, performance and versatility for both research and OEM applications. With peak UV QE at 32% and peak VIS-NIR at 58%, Sincerity™ offers a relatively flat response and addresses multiple applications. In the near-NIR, this detector is a lower cost alternative to a deep depleted CCD, with no etaloning. The dynamic range of the Sincerity™ camera is unmatched at 92.5 dB (in image mode) and the CCD's tall format is an excellent fit for imaging spectroscopy applications. Its flexible design allows our OEM-dedicated team to quickly adapt the camera for industrial requirements, ranging from alternate CCD chips to electronics customizations.

Feature	Spectroscopy Benefits for Research and OEM
Unmatched Readout Noise	Unmatched Readout Noise at a speed of 45 kHz
Unparalleled Linearity & Dynamic Range	Measured values on each camera exhibit outstanding performance: 92.5 dB typical dynamic range with <0.4% non linearity.
Deep Thermoelectric Cooling (Minimum)	–60°C @ +25°C ambient yields very low dark current
Ultra-Compact size	Ideal for use on microscopes and OEM integration
Lifetime Vacuum Warranty	All metal sealed technology provides a maintenance-free permanent vacuum.
PC Interface	USB 2.0 high speed with 100% data integrity. No controller box.
Ruggedized Connectors	Maintains overall system integrity in industrial environments.
Open Electrode CCD Technology	Broad spectral coverage with QE of 27% at 250 nm and 55% at 800 nm.
Scientific Grade 1 CCD with 6.7 mm Height	Ideally suited for low light level detection in a variety of spectroscopic applications. We offer other CCD formats and sensor types for OEM volumes.
Flexible Input & Output Trigger Interface	Experiment synchronization with External Trigger In & TTL Shutter Out with programmable edge triggering
LabView VIs and SDK available.	Flexible software to integrate a Sincerity CCD into existing apparatus or as an OEM component. Contact us for a Linux driver.



# Specifications

## SYNCERITY™ 1024 × 256 – Open-Electrode


<b>CCD Sensor Format</b>	1024 × 256
<b>Quantum Efficiency at 20 °C (See QE curve below)</b>	27% at 250 nm 31% at 300 nm 42% at 550 nm 58% at 750 nm 55% at 800 nm 12% at 1000 nm
<b>Pixel Size</b>	26 μm × 26 μm
<b>Image Area</b>	26.6 mm × 6.7 mm, 100% fill factor
<b>Deep Thermoelectric Cooling</b>	–60 °C @ +25 °C ambient or –50°C @ +40 °C ambient Yields low dark current suitable for most OEM and some Research applications
<b>Single Pixel Well Capacity</b>	200 000 e <sup>-</sup> /pixel (Minimum)
<b>Serial Register Full Well Capacity</b>	1 000 000 e <sup>-</sup> /pixel (Typical Output Register Saturation)
<b>Scan Rates</b>	45 kHz and 1 MHz
<b>Readout Noise (at 45 kHz and at –60 °C) *1</b>	4.7 e <sup>-</sup> (Typical) to 7 e <sup>-</sup> (Maximum)
<b>Readout Noise (at 1 MHz and at –60 °C) *1</b>	17 e <sup>-</sup> (Typical) to 20 e <sup>-</sup> (Maximum)
<b>Maximum Spectral Rate</b>	27 Hz at 45 kHz scan rate 278 Hz at 1 MHz scan rate
<b>Digitization</b>	16 bit ADC
<b>Dynamic Range (Typical for Single Pixel) *2</b>	42 550:1 (92.5 dB providing >15 bit effective dynamic range)
<b>Non Linearity (Measured on Each Camera)</b>	< 0.4% at 45 kHz – Linearity better than 99.6% < 0.8% at 1 MHz – Linearity better than 99.2%
<b>Dark Current at –60 °C *3 (Note that pixel size = 26 μm)</b>	0.0052 e <sup>-</sup> /pixel/sec (Typical) equivalent to 0.0020 e <sup>-</sup> /pixel/sec for a 16 μm pixel size equivalent to 0.0031 e <sup>-</sup> /pixel/sec for a 20 μm pixel size
<b>Software-Adjustable Gains</b>	1–12 e <sup>-</sup> /count
<b>Environmental Conditions</b>	<ul style="list-style-type: none"> <li>o Operating Temperature 0°C to 40°C ambient</li> <li>o Relative Humidity &lt; 70% (non-condensing)</li> <li>o Storage Temperature –25°C to 50°C</li> </ul>
<b>Weight</b>	1.769 kg (3.90 lb)
<b>Dimensions</b>	Refer to mechanical drawings herein
<b>Power Requirements</b>	
<b>AC-DC Power Supply (Provided)</b>	90–264 VAC, 47–63 Hz
<b>Recommendation for OEM Supplying Camera Power Directly:</b>	<ul style="list-style-type: none"> <li>• Pin: +9 V, ± 5%, 6.44 A maximum</li> <li>• Regulation: +8.55 V<sub>min</sub>, +9 V<sub>typ</sub>, +9.45 V<sub>max</sub></li> <li>• Ripple &amp; Noise: 200 mV<sub>pp</sub> maximum</li> </ul>
<b>Minimum Computer Requirements:</b>	<ul style="list-style-type: none"> <li>• 3.0 GHz single core or 2.4 GHz multi-core processor</li> <li>• 2 GB RAM</li> <li>• 32 bit or 64 bit compatible</li> <li>• 500 MB free hard disk space (additional disk space may be required depending on data storage needs)</li> <li>• USB 2.0 High Speed Host Controller capable of sustained rate of 40 MB/s</li> <li>• Windows (XP, Vista and 7)</li> </ul>

\* All specifications subject to change without notice.

### Footnotes:


1. Entire system noise measured for a single pixel
2. Dynamic range is defined as: Full Well / Readout Noise and is measured at 45 kHz
3. Averaged over CCD area, but excluding any regions of blemishes.

# Quantum Efficiency




# Connecting to Syncerity

<b>Power Interface:</b>	
Connector Type:	PDP-40, Mini PWR DIN, 4-Position, STR Plug
<b>Camera Interface:</b>	
Connector Type:	USB standard Type B
<b>Sync I/O Connectors:</b>	
Connector Type:	SMB
Input Jack:	TTL IN (EXT TRIG IN)
Output Plug:	TTL OUT (SHUTTER OUT)



# Product Drawings

Dimensions in inches (mm)



## Ordering Information

### Syncer-1024x256-OE Syncerity TE Cooled CCD Camera Includes:

- USB 2.0 Camera Head
- Power Supply AC-DC
- USB Cable
- CD Manual

### Optional:

- Shutter Driver (SDrive-500 Shutter Control Unit with cable)
- CCD Shutter
- TTL (IN) Trigger Cable
- Printed Manual

# Accessory Information:

## SDrive-500 Shutter Control Unit

The SDrive-500 Shutter Control Unit is designed to control the activation/opening of an electro-mechanical shutter during the interval when the CCD of a HORIBA Scientific detection system is being exposed to light.

Specifications	
System Parameter	Units/Description
Environmental	
Environmental Conditions	Storage Temperature from $-25^{\circ}\text{C}$ to $+50^{\circ}\text{C}$ Operating Ambient Temperature Range $+25^{\circ}\text{C} \pm 5^{\circ}\text{C}$ Relative Humidity $\leq 70\%$ non-condensing
TTL Level Electrical Input Signals	
Trigger Input #1	TTL level input signal triggering shutter activation on the rising edge.
Trigger Input #2	TTL level input signal triggering shutter activation on the rising edge.
Shutter Override	TTL level switch located on the unit's front panel providing a manual override activation (i.e., opening) of the shutter.
Shutter Output Excitation Drive	
Shutter Coil Resistance	$12\ \Omega$
Shutter Pulsed Voltage to Open	+48 VDC
Shutter Hold Voltage	+5 VDC
Operating Frequency	40 Hz maximum rep rate
Input Power Requirements	
Input Line Voltage	85–264 VAC continuous / universal
Input Line Frequency	47–63 Hz
Input Power	55 W max
Mechanical Dimensions	
Dimensions (L×W×H)	7.24 in × 4.50 in × 3.15 in
Weight	1.3 kg (2.9 lbs)

• All specifications subject to change without notice.

