



Large Capacity Laboratory Refrigerators

MPR-1410 MPR-1410R MPR-721 MPR-721R



Stable and reliable refrigerated environment for exacting laboratory requirements. Stable temperature environment for clinical research, pharmaceutical and industrial use. Adjustable shelves and wide range setpoint for varying laboratory applications including chromatography and pharmaceutical storage.

> SANYO is committed to developing green technologies that provide energy efficiency resulting in lower operational costs with less impact on the environment.

Model: MPR-1410



Large Capacity Laboratory Refrigerators: MPR Series

Temperature Stability

SANYO's temperature control system, with thermistor temperature monitor and microprocessor control reliabably maintains cabinet temperature at the set level and is unaffected by outside ambient temperature. Forced air circulation ensures that the cabinet temperature returns to the set point quickly after door openings and remains uniform throughout the cabinet.

Wide Temperature Range

With a temperature range of 2°C to 23°C, SANYO's MPR Series refrigerators are ideally suited for many tests that require a stable, cool temperature or general cold storage.



Easy Access Front Control Panel



Speedy and Powerful Refrigeration

For frequent door openings, SANYO's MPR Series refrigerators are equipped with powerful, hemetically sealed compressors. These purpose-built compressors ensure superior pull-down characteristics and precise temperature control.

Easy-to-Manage Layout

The interior layout flexibility of the MPR Series refrigerators makes them ideal for running experiments that require stable, cool conditions, as well as general cold storage.

Standard Alarm and Safetly Features

SANYO's MPR Series refrigerators are fitted with alarms and flashing indicator LED to warn of high and low temperature conditions. In the event of an irregular rise in cabinet temperature, the heater automatically shuts off and forced air circulation brings the temperature down. Door locks are standard to safegaurd valuable contents.

HFC Refrigerant & CFC-Free Insulation

SANYO biomedical equipment is designed for low environmental impact. The MPR-721 and 1410 series use HFC refrigerant R-134A, and foamed-in-place insulation is CFC free.



Large Capacity

With a modular width of 30.3 inches (770mm) the MPR-721 offers capacity of 24.1 cu.ft (684 liters) [MPR-721R: 23.7 cu.ft. (671 liters)], while the 56.7 inches (1440 mm)-wide MPR-1410 offers capacity of 48.4 cu.ft. (1370 liters) [MPR-1410R: 48.2 cu.ft. (1365 liters)]. The interior is spacious enough to accomodate column chromatography apparatus or large volumes of reagents, test samples and biologicals.

Large Fans

The 4.7"-diameter fan ensures an even temperature throughout the cabinet (MPR-1410R models have a double flow system with two fans). Heat spots from powered test apparatus are minimized and pull-up characteristics after door openings are outstanding.



Stable & Reliable Refrigerator Environment

MPR-1410 Performance Data



Adjustable Shelves (MPR-721 & MPR-1410)

The shelves can be arranged to accommodate tall apparatus such as fraction collectors. These shelves are deep enough [24.4 inches (620mm) front to back] and strong enough (110lbs load for the MPR-721, 88lbs for the MPR-1410) to hold most apparatus.



Drawer Type (MPR-721R & MPR-1410R)

The "R" models are fitted with pull-out drawers. With a profile of 3.9 inches (100mm) and 20.8 inches (530mm) front to back, these drawers are deep enough to hold large bottles or reagent kits. They also allow convenient, space-efficient storage and management of patient medications and other items. * The MPR-721R is shipped with five drawers, while the MPR-1410R is shipped with ten drawers, five in each half of the cabinet.

Features

It is common knowledge that product should be stored in stable conditions below ambient temperature. Domestic refrigerators are capable of storage at 4°C, but they suffer from the following drawbacks:

- 1. Temperature varies each time the door is opened.
- 2. Temperature rises during defrost cycle.

Cabinet temperature is easily affected by ambient temperature with the risk of contents freezing if ambient drops below 0°C.
Temperature setting by dial is inaccurate (no digital temperature indication).

SANYO has incorporated solutions for all these problems into its biomedical refrigerators, which have been well received in hospitals and research facilities around the world.

MAIN FEATURES

1. Stable, uniform and controlled cabinet temperature is unaffected by outside temperature.

2. Cycle defrost allows defrosting without rises in cabinet temperature.

3. Standard alarm & safety features prevent irregular temperature flucuations in cabinet.

User Friendly Design & Double Door Function

Easy to read and operate control panel features a full array of alarm and safety functions. Through wall access ports are standard to allow access for power cords. Aesthetically attractive exterior color blends in most lab color schemes. The "catch-free" rounded corners are safe and attractive. Doors open smoothly and close automatically. Some laboratory



refrigerators have solid doors for temperature stability, while others have glass doors for easy viewing. SANYO MPR Series refrigerators give you the best of both versions. The broad, solid frames with effective gaskets ensure excellent temperature stability, and the double-pane glass windows offer an excellent view of the interior without compromising temperature stability. Optional panels are available to block out light.

Cycle Defrost

In the 4°C range, frost buildup on the evaporator is inevitable. This can affect the performance of the heat exchanger. SANYO has solved this problem with a cycle defrost and evaporator temperature sensor system. This system runs automatically so there is no need to turn off power for defrosting. Temperature rise during defrost is minimal. The evaporation heater also doubles as protection against drops in cabinet temperature caused by a low ambient temperature.



MPR-721 Performance Data



MPR Series

Model	Shelves Drawers
MPR-1410	8 Wire Shelves
MPR-1410R	10 Drawers
MPR-721	4 Wire Shelves
MPR-721R	5 Drawers

New concept features have been added to the MPR-721 and MPR-721R. A filterless condenser cuts down on cleaning maintenance. Also, all glass doors now include a protective film to prevent shattering in the event of accidental collisions.



MPR-1410 Effective capacity 48.4 cu.ft.



MPR-1410R Effective capacity 48.2 cu.ft.



MPR-721 Effective capacity 24.2 cu.ft.



MPR-721R Effective capacity 23.7 cu.ft.

Dimensional Data









MPR-1410



MPR-1410R





MPR-721



MPR-721R

Features	Trigger	Type of alarm/response			
High temperature alarm	1410R: 5°C or more above setpoint 721R: 2°C to 14°C (selectable) above setpoint	Flashing LED & audible alarm after 15 minutes			
Low temperature alarm	1410R: 5°C more below setpoint; 0°C or lower 721R: 2°C to 14C (selectable) below setpoint; 0°C or lower	Flashing LED & audible alarm			
Over heating protection	Interior temperature rises to 40°C	Fan motor & heater OFF			
Temperature lock	Lock ON	Set temperature cannot be altered			
Memory backup	Power failure	Memory of settings			
Door ajar	Door not closed	Door ajar lamp ON			
Self diagnostics	Sensor open or short	Error code shown (E1, E2)			
Remote terminal	Operation can be checked from office or control room away from the lab				
Temperature recorder	Optional, Model # MTR-C954				

Alarm & Safety Features

Specifications



MPR-1410, 1410R, 721 and 721R

Model N	0.	MPR-1410	MPR-1410R	MPR-721	MPR-721R		
Exterior	dimensions	56.7" x 36.2" x 76.8"		30.3" × 36.2" × 77.0"			
(W x D x	H)	1440 x 830 x 1950mm		770 x 830 x 1955mm			
Interior c	limensions	52.0" x 28.0" x 59.1"		25.6" x 28.0" x 59.1"			
	п) 			650 x /10 x 1500mm			
Effective	сараситу	48.4 cu.π. (1370 liters)	48.2 cu.π. (1365 liters)	24.2 cu.π. (684 liters)	23.7 cu.π. (671 liters)		
Exterior		Acrylic finish baked on zinc galvanized steel					
Interior			Acrylic finish baked or	n zinc galvanized steel			
Doors		2 x double pane glass doors, self closing		1 x double pane glass doors, self closing			
Insulatio	n	CFC-tree rigid toamed-in-place polyurethane					
Shelves/	Drawers	8 x polyethylene-coated wire shelves	10 x coated steel drawers.	4 x polyethylene-coated wire shelves	5 x coated steel drawers.		
		Max. load: 40kg/shelf	Handles with card holder.	Max load: E0kg/shalf	Handles with card holder.		
			Max. load: 40kg/drawer	IVIAX. IUAU. SUKY/SHEII	Max. load: 40kg/drawer		
Access p	orts	3 x 30mmø (2 in sides, 1 in cabinet top)					
Locks	2 x cylinder type		1 x cylinder type				
Casters		4					
Cooling	ig method Forced air circulation with double fan		Forced air circulation with single fan				
Compres	ompressor 300W hermetic rotary		220W hermetic rotary				
Evaporat	tor	Fin tube type					
Condens	Fin tube type		Wire tube type, Filterless design				
Refrigera	int	R-134a (HFC)					
Defrost r	nethod	Forced type (cycle defrost system), fully automatic hot pipe for automatical evaporation for drain water					
Defrost h	neater	257W		153W			
Power	Voltage	11	115		115		
source	Hz	60		60			
	Phase	1ø		1ø			
	Amps	3.7		2.4			
Max. hea	at release	1565btu		1440btu			
Tempera	ture control	Microprocessor controlled on/off type. Thermistor sensor Temp. control range: 2°C to 14°C (ambient temp. –5°C~0°C no-load) 2°C to 23°C (ambient temp. 0°C~35°C no-load) (1°C step)					
Tempera	ture display	Digital type (1C step)					
Safety fe	atures	High temp. protection circuit, low temp. protection circuit, temperature lock, self diagnostics, memory backup (nonvolatile memory).					
High terr	ıp. alarm	When temp. 5°C or more above set level, flashing LED & buzzer after 15-minute delay Flashing lamp when temp. 5°C or more above set point, buzzer after 15-minute delay					
Low tem	p. alarm	When interior temp. drops to 0C or lower, flashing LED & buzzer					
Door aja	r alarm	Flashing lamp, buzzer after 2-minute delay					
Interior la	amp	40W x 1, fluorescent lamp 20W x 1, fluorescent lamp			rescent lamp		
Net weig	jht	509.3 lbs (231kg)	588.63 lbs (267kg)	368.2 lbs (167kg)	410 lbs (186kg)		

*Appearance and specifications are subject to change without notice.

SANYO offers a wide range of high quality validation services for all of our equipment. These services include on site validation, custom validation, custom validation support packages, factory acceptance testing, critical utility qualifications, NIST traceable calibration. Choosing SANYO as an equipment supplier and validation consultant can greatly reduce the time and cost involved with setting up and operating new equipment.



SANYO North America Corporation Biomedical and Environmental Solutions 1300 Michael Drive, Wood Dale, IL 60191 USA Toll Free USA 800-958-8442 • Fax 630-238-0074 www.sanyobiomedical.com

SANYO Canada, Inc. 201 Creditview Road, Woodbridge, Ontario L4L 9T1 905-265-4100 • Fax 905-265-4101

LR-MPR/LabsSeries-7.09