CERTOMAT® Tplus

Inkubations-Schüttelschrank

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

CERTOMAT Tplus: ca. 200kg

Untergestell Tplus: ca. kg

The reproduction, distribution and utilization of this document as well as the
 creation of new documents by means of reprographic or electronic means is
 prohibited without the written permission of Sartorius. All rights reserved in the
 extent of the grant of a patent, utility model or design.
 Änderungen nur über CAD!

Index		Field		Changes/Änderung	
HINWEIS: RoHS-Konformität beachten ! NOTICE: please note RoHS-conformity !					
General tolerances/ Allgemeintoleranz		Version	Changes Änd.	Date Datum	Name
ISO 2768 - m H		-	-	-	-
Kabelängen (cable lengths) ±10%		-	-	-	-
<input type="checkbox"/> Maße nach DIN 121 überprüfbar <input type="checkbox"/> measurements to DIN 121		Drawn by	28.02.07	Strietzel	Part Name/Benennung
<input type="checkbox"/> Holes according to <input type="checkbox"/> unbracketed workpiece edges according to ISO 13715		Checked by	28.02.07	Sandri/Strie	Certomat T plus Aufstellplan
		Approved by	-	-	
		sartorius MT-Göttingen			Drawing No./Zeichn.-Nr. B370380000101 - A1
Replaces Ers. d.		Replaced by Ers. d.			Scale/Maßstab . . . Sheet Blatt . . .

CERTOMAT® Tplus

Incubation Shaking Cabinet

Operating Manual

Introductory Notes

This Operating Manual for the Incubator Shaking Cabinet CERTOMAT® Tplus refers to the design and equipment features of this device at the date of print given with the Rev.No.. The information included in this Operating Manual has been carefully estimated and checked. However, no guarantee can be given with respect to completeness and correctness of this information. Sartorius Stedim Biotech GmbH reserves the right to modify the equipment and to change this operating manual without notice.

Due to the continuous development of the units it may be possible that different building components or their use are not described. There may be information which does not apply to your device. If you find missing, misleading or wrong information or if you should wish additional information about different components please send us a note.

The CERTOMAT® Tplus is just an example of laboratory equipment products from Sartorius Stedim Biotech GmbH. Our range of equipment includes shakers and incubator shakers, homogenizers, centrifuges, ultrapure water systems, air samplers and fermenters/bioreactors.

You will find further information on our homepage www.sartorius-stedim.com.

Upon request, we will be pleased to inform you about our complete product range.

Sartorius Stedim Biotech GmbH
August-Spindler-Straße 11
37079 Göttingen, Germany
Telephone +49.551.308.3118
Telefax +49.551.308.3918
info@sartorius-stedim.com
www.sartorius-stedim.com

EG-Statement of Conformity

With the Statement of Conformity, which is attached or which accompanies the equipment, Sartorius Stedim Biotech GmbH confirms the conformance of the CERTOMAT® Tplus with the specified guidelines.

Intended Use

Incubation shaking cabinets are power-driven machines (acc. to BGR 500 Chap 2.11 Part 3) in which liquids or nutrient solutions are shaken with cultures of biological origin and which are kept at a selected incubation temperature and, hence, are also only intended for this purpose. Any other or further use is considered to be unauthorized. Sartorius Stedim Biotech GmbH does not accept liability for any damage resulting from unauthorized use.

Authorized use also includes:

- paying attention to all instructions from the operating manual and
- adherence to service, cleaning and maintenance regulations.

General Safety Notes

- Laboratory seniors or the persons responsible for the use of the equipment must inform the personnel about potential risks of operating the CERTOMAT® Tplus. Furthermore they have to guarantee that the relevant safety instructions are observed. The personnel must have the necessary qualification for the operation of the equipment and its handling.
- Please take special care regarding the operating information. Improper use of the unit can thus be avoided and you retain full claim for guarantee. Please contact us if you have any questions.

Steps to be carried out with special care or special aspects or special issues to be considered are marked as this paragraph.

Important safety notes in this manual are marked with this symbol. Such instructions apply in addition to the other safety instructions for the field of application and for the working place.

- Do not use devices with the wrong power rating.
- Disregard of the safety instructions may for example result in damage to the equipment or may cause other material and eventually also person damage.

Guide through this Manual

This Operating Manual is organized systematically. Depending on your knowledge of the device you can start from the beginning or go directly to the chapter of your interest.

Chapter 1 provides an overview of the construction and design of the CERTOMAT® Tplus.

Chapter 2 provides information on the requirements regarding the placement of the device. In addition, it describes the installation and connection as well as equipment configuration, loading and mounting of trays and universal assembly systems.

Chapters 3 to 11 give information on the commissioning and operation of the device as well as information on measures for operational malfunctions. If you are already familiar with the device and it has been installed and is ready for operation, you can go directly to this part of the manual.

Chapter 12 contains technical data and ordering information. This chapter includes information on optional equipment and the range of accessories for the shaker in particular.

Contents

Introductory Note

EG Statement of Conformity

Intended Use

General Safety Notes

Guide through this Manual

1	Design and Function	5	5	"Main" Function	22
1.1	Design and Construction Features	5	5.1	Functional Elements of the "Main" Menu	22
1.2	Available Equipment and Accessories	5	5.2	Starting an Incubation from the "Main" Menu	22
1.2.1	Shaking Drives	5	5.2.1	Entering the Temperature	22
1.2.2	Trays	5	5.2.2	Entering the Speed	23
1.2.3	Illumination Unit	5	5.2.3	Switching the Working Light on or off	23
1.2.4	Grid for Petri Dishes	5	5.2.4	Switching the Photo Illumination on or off	23
			5.2.5	Switching the Cooling on or off	23
2	Delivery and Installation	6	5.3	Editing individual Values of an Incubation from the "Main" Menu	24
2.1	Delivery and Checking of Completeness	6	5.4	Ending an incubation from the "Main" menu	24
2.2	Installation and Laboratory Connection	7	5.4.1	Switching off the heating	24
2.2.1	Placement at the Working Place	7	5.4.2	Switching off the Shaker (Speed 0 rpm)	24
2.2.2	Power Supply	7	5.4.3	Switching the Working Light, Photo Illumination and Cooling on or off	24
2.3	Configuration of the CERTOMAT® Tplus and Optional Equipment	8	6	Main "Controller" Function	25
2.3.1	Shaking Drive with Tray Support	8	6.1	Functional Elements of the Main "Controller" Menu	25
2.3.2	Illumination Unit (option)	8	6.2	Access to the Controller	25
2.3.3	Connection of the Optional Cooling Devices	8	6.2.1	Entering the Set Temperature Value	25
2.3.4	Installation of the Grid for Petri Dishes (optionally)	8	6.2.2	Setting the Alarm Limits	26
2.3.4.1	Installation of the Mounting Rails (4 pieces)	9	6.2.3	Select "Controller Operating Mode".	26
2.3.4.2	Installation of the Angle Bracket	9	7	Main "Batches" Function	27
2.4	Equipment, Mounting and Loading of a Tray	10	7.1	Functional Elements of the Main "Batches" Menu	27
2.4.1	Equipping the Tray	10	7.2	Example of Creating a Program	27
2.4.2	Adjustment of the Counter Weight to the Intended Load	11	7.3	Entering Program Data	29
2.4.3	Mounting and Dismounting a Tray	12	7.4	Specifying the Start Time of a Program	33
2.5	Mounting of Trays and Mounting Sets	13	7.4.1	Starting the Program immediately	33
2.5.1	Assembly and Loading of Finally Equipped Sets of Trays	13	7.4.2	Starting the Program at a specified Time	34
2.5.2	Mounting of Test Tube Racks	13	7.5	Stopping / Cancelling a Program	37
2.5.3	Mounting of the Universal Mounting System	14	7.6	Editing the Program	38
3	Operating Information	15	8	Main "Settings" Function	40
3.1	Safety Notes	15	8.1	General Information	40
3.2	Equipment, Mounting and Loading of a Tray	16	8.2	Functional Elements of the Main "Settings" Menu	40
3.2.1	Equipping the Tray	16	8.3	Operating Menus and Input Window	41
3.2.2	System behaviour at the Start	17	8.3.1	"System Parameters"	41
4	Basics of the Operation of the "Touch Display"	18	8.3.2	"PV Ranges"	41
4.1	Layout of the main screen	18	8.3.3	"Calib. Temperature"	41
4.1.1	Header	18	8.3.4	"Service"	41
4.1.2	Working Area	19	8.3.5	„Manual Operation“	41
4.1.3	Footer	19			
4.2	Display of the Functional Elements	20			
4.3	Overview of the Main Function Keys	21			

9	"Remote Function"	42
9.1	Functional Elements of the "Remote" Function	42
10	Main "Alarms" Function	42
10.1	Functional Elements of the Main "Alarms" Menu	42
10.2	Confirming and/or deleting Alarms	43
10.3	Overview of Process Alarms	43
10.4	Process Messages	43
10.5	Abbreviations	44
10.5.1	Process Values	44
10.5.2	Digital Inputs	44
10.5.3	Digital Outputs	44
11	Appendix	45
11.1	Operation of the Internal Cooling Unit	45
11.2	Connection of External Equipment	45
11.2.1	Analog Exits (ANALOG OUT)	45
11.2.2	Collective Alarm	45
11.3	Maintenance and Trouble-Shooting	46
11.3.1	Electrical Malfunctions	46
11.3.2	Malfunctions of the Shaking Drive	46
11.3.3	Malfunctions of the Optional Cooling Unit	46
11.3.4	Handling and Eliminating Faults in the DCU System	46
11.4	Maintenance	47
11.4.1	Maintenance of the Illumination Unit	47
11.5	Cleaning	47
11.6	Warranty Regulations and Service	48
11.7	Password System	48
11.8	Installation Plan, Interface Configuration	48
12	Technical Data, Ordering Information	49
12.1	Technical Data	49
12.1.1	Mechanical Design	49
12.1.2	Mains Supply	49
12.1.3	Operating Data	49
12.1.4	Analog Out	49
12.1.5	Collective Alarm	49
12.1.6	Option Cooling Unit	49
12.2	Ordering Information	50
12.2.1	Configuration of the Shaking Cabinets	50
12.2.2	Optional Equipment	50
12.2.3	Accessories	51
12.3	Safety Instructions Laboratory Shakers and Incubation Shaking Cabinets	57
12.4	EG Declaration of Conformity	57
12.5	Declaration of Decontamination	57
12.6	Information and Instructions on Disposal and Repairs	57
12.7	Dimension sheet CERTOMAT® Tplus	57
12.8	Prospectus CERTOMAT® Tplus	57

1 Design and Function

1.1 Design and Construction Features

The CERTOMAT® Tplus is an easy to operate, efficient and durable incubator shaking cabinet for universal use in the various tasks of biological and chemical laboratories, for example.

The basic equipment of the CERTOMAT® Tplus includes one tray support for shaking trays, the drive unit, and the heater. The drive consists of a motor of external rotor type with brushless drive and a noiseless, robust Poly-V belt-drive. Shaking drives with different orbits are available. The orbit can be modified by our service on site. The different shaking speeds can be adjusted continuously. The shaker can be equipped with trays of different size and numerous clamps, holders as well as mounting systems of the CERTOMAT® - shaker program.

Two large doors facilitate the fixing and removing of the shaking containers or the changing of the tray and they allow the visual control. For incubations in the temperature range near to and below room temperature a cooling unit can be installed. More information concerning the equipment configurations you will find in chapter 12 Technical Data, Ordering Information. With this information you can specify the optional equipment of the shaking cabinet according to your specific demands.

The incubation chamber of the cabinet is made of stainless steel. At the bottom a spill tray made of stainless steel is mounted. If containers get damaged during shaking operation this prevents spilled media from escaping unchecked. The spill tray can be easily removed from the incubation chamber and then cleaned. This design also facilitates the efficient cleaning of the shaking cabinet.

You can use the shaking cabinet on the floor or on the table and you can stack 2 – 3 shaking cabinets. For stacking 2 cabinets a support frame is available. That means optimal arrangement for operation.

1.2 Available Equipment and Accessories

1.2.1 Shaking Drives

The operating terminal and the control system are included in the panel on the right side of the cabinet. The shaking drive is mounted below the spill tray inside the cabinet. The drive is similar to the drive used for the bench top shaker CERTOMAT® RM. The drive specifications are as follows:

- Orbital motion Ø 25 mm and 50 mm (depending on the version delivered; can be changed to alternative amplitude on site by authorized service representatives of Sartorius AG).
- Shaking speeds: 40 ... 400 ¹/min
- Rotation speed accuracy: max. ± 1% of final value
- The drive has a patented, variable mass compensation.

1.2.2 Trays

Trays can be delivered with sets of fixed clamps for Erlenmeyer flasks of different sizes (type E, F), or they are available for individual mounting as universal trays type EU or FU.

- Dimensions:
 - type E/EU: W x D about 420 x 420 mm ¹
 - tray type F/FU: W x D about 800 x 420 mm ¹

The universal trays can be equipped with clamps for Erlenmeyer flasks of different sizes, hinged test tube racks and universal mounting system. The universal mountings allow to use various vessels, such as beakers, bottles, or separation funnels. In addition, non slip mats and adhesive tape are available.

1.2.3 Illumination Unit

- The Tplus can be ordered with an optional illumination unit comprising 5 lamps of 18 W each. Each lamp can be switched off separately. The mounting of the illumination unit is described in chapter 2.3.2 Illumination Unit (option).

1.2.4 Grid for Petri Dishes

- For static incubations, e. g. in Petri dishes, a grid for petri dishes is optionally available. The mounting of the tray insert is described in chapter 2.3.4 Installation of the Grid for Petri Dishes (optionally).

¹) Trays of modified design for actual locking; rebuilding of old trays for this locking on request. Please contact the Sartorius customer service.

2 Delivery and Installation

2.1 Delivery and Checking of Completeness

The configuration and the extent of delivery of the incubation shaking cabinet depend on the individual order. The shaking drive, and if ordered, the option „cooling unit“ will be preinstalled in the factory prior to delivery according to the customer’s demands. Each CERTOMAT® Tplus is delivered after a thorough functional test. For transport to the working place, and prior to installation and start-up of the unit you should check the following:

1. A suitable transport route through the building and appropriate transportation means (lifting trucks, etc.) should be available. The floor, the doors and elevators must be large enough for safe transport.
2. The working place must be suited for the shaking cabinet. It must offer sufficient space for convenient placement and it must compensate the weight of the completely mounted shaking cabinet:
 - Single cabinets or a stack of two cabinets can be placed on a support frame. This will ease operation of the controls.
 - Dimensions of the support frame:
1150 x 200 x 710 mm (W x H x D)
 - Weight of CERTOMAT® Tplus:
approx. 230 kg including load

For devices with integrated cooling unit (UHK versions), the following applies for incubations below the ambient temperature:

If the machine is operated continuously for a long period in conditions of relatively high humidity, ice can form on the vaporiser.

The formation of ice on the vaporiser reduces the air circulation.

This can activate the overheating safety cut-out and switch off the machine.

Allow the machine to cool down before turning it on again.

- For placement of one or 2 cabinets in stack arrangement we recommend to place the (lower) cabinet onto the optional support frame. This makes the operation easier.
 - Up to three cabinets can be stacked. In this case the working place must be steady enough to carry the weight of all cabinets including all accessories and the load. Also the dimensions of the devices must be taken into consideration.
3. The working place must be able to compensate vibrations which are caused by the shaking operation. Make sure that surrounding equipment cannot be disturbed by the Tplus. Please observe the notes for the limitation of the shaking speeds in chapter 3.

The action of the loaded trays causes vibrations, which will change with load and shaking speed. These vibrations cannot be fully compensated by the adjustable counter-balance weight. Especially at heavy loads (>10 kg) and high shaking speeds resonance effects may occur and vibrations may be transmitted to the working place.

4. The power supply of the laboratory must correspond to the technical specifications of your CERTOMAT® Tplus -see the labels on the unit. Compare the existing line voltage with the data on the nameplate of the device.
5. The shaking cabinet has been designed for normal laboratory conditions, see information about the environmental conditions in the supplement. Under special environmental conditions or when using aggressive media you have to test the suitability of all parts. Malfunction and damage, such as corrosion by aggressive media, are not subject of our guarantee.

If the working place is not suitable for the CERTOMAT® Tplus and the unit does not correspond with your power rating, you must not start and operate it.

6. Please check the condition of your device and the accessories provided. All components specified in the order are premounted or are enclosed in the delivery. The delivery should be complete and all parts undamaged.
- If the CERTOMAT® Tplus has the wrong mains specifications, in case of malfunction of the unit or if any parts are damaged or are missing inform your Sartorius Stedim Biotech GmbH representative as soon as possible.

2.2 Installation and Laboratory Connection

2.2.1 Placement at the Working Place

You can stack the CERTOMAT® Tplus. Up to three devices can be placed on top of each other. As far as the information below refers to a shaking cabinet they apply also to your unit.

1. Place the CERTOMAT® Tplus in such a way, that it is easily accessible and does **not** interfere with other devices. The right part of the machine should be positioned in an approx. distance of 50 cm from other devices or from the wall. So the side door is accessible and you will not have to reposition the CERTOMAT® Tplus, if you want, for instance, to open the side door. A minimum wall distance of 15 cm is necessary for sufficient ventilation.
If the CERTOMAT® Tplus includes the optional cooling unit (UHK version), a clearance from the wall of at least 20 cm is required.
- For placement of one or 2 cabinets in stack arrangement we recommend to place the (lower) cabinet onto the optional support frame. This makes the operation easier. Each unit can be operated separately.
2. To stack to 2 to 3 cabinets
 - remove the plastic caps at the upper cover of the cabinet, which is to be placed below.
 - replace the feet of the intended upper cabinet by the locating pins included in the delivery.
 - if necessary, the eyelets provided can be screwed into the upper side of the top cabinet so that the cabinet can be easily lifted.

When screwing in the lifting eyebolts use the plastic washers provided. The washers prevent the cover plate from being scratched.

3. Place the cabinets carefully on top of each other. Take care of the weight of the cabinets. Only use suited lifting tools or devices. The locating pins should easily fit into their borings in the upper cover of the lower cabinet.
4. Carefully level the incubation shaker or the stacked units.

Unbalanced placement of the cabinet can cause uneven rotation of the tray support. The drive can be subject to quick wear and tear or be damaged, respectively.

Especially at heavy loads of the shakers and with high shaking speeds, the stack arrangement of 2 ... 3 cabinets can start to wobble during operation.

- Carrying handles are provided for easier handling of the device.

Handling:

1. Switch off the CERTOMAT® Tplus and disconnect it from the power supply.
2. Open the side door on the right and push the assembled handle into the two drill holes in the frame. Secure the carrying handle on both sides against slipping with the enclosed split pins.
3. Push the left side wall slightly upwards so that it slides out of the holding rails. Remove the side wall and put it to one side.
4. Push the assembled handle into the two drill holes in the frame. Secure the carrying handle on both sides against slipping with the enclosed split pins.

The CERTOMAT® Tplus can now be moved.

CAUTION:

The device is very heavy! It can only be moved by at least 4 persons.

5. When the work is ended, pull out the split pins and remove the handles.
6. Close the right side door and put the left side door onto its hinges again.
7. Only after the device is fully closed once more can it be connected to the supply voltage again.

2.2.2 Power Supply

- You will require or need to prepare a laboratory grounded socket offering 230 V, 50 Hz or 115 V, 60 Hz, depending on the country standard. Check on the type plate on the cabinet whether it has the correct power ratings.
Never connect the CERTOMAT® Tplus to power supplies of wrong ratings.
- If the plug at the power supply cable does not fit to your sockets, the service technician of Sartorius AG can attach a suitable plug.

2.3 Configuration of the CERTOMAT® Tplus and Optional Equipment

2.3.1 Shaking Drive with Tray Support

- The shaking drive will be mounted in the factory prior to delivery in the shaking cabinet and will be ready for use after assembly of the tray and loading with shaking containers. Rebuilding to a different orbit must only be done by authorized and qualified service personnel. Please contact your service representative or the service department of Sartorius AG, if necessary.

2.3.2 Illumination Unit (option)

- Push the illumination unit in the prepared guiding bars. You have to release the 4 fastening screws of the pull bars first. Push the illumination unit into the device until the cross hole in the pull bars becomes visible. Insert the delivered screws in the cross holes and protect the illumination unit in this way from sliding out. Fasten the 4 fastening screws carefully again. Plug the plug for the mains supply into the prepared socket in the incubation chamber (refer to fig. 1). Secure the plug (refer fig. 2) by fastening the connecting nut in the socket (refer fig. 3). The illumination unit is then ready for operation. If required, individual lamps can be manually switched off to reduce the light intensity (refer fig. 4).

Fig. 1: Socket for the optional illumination unit

Fig.. 2: Plug of the optional lighting unit

Fig. 3: Affiliated illumination unit

Fig. 4: Switch for fluorescent tubes

2.3.3 Connection of the Optional Cooling Devices

- If your CERTOMAT® Tplus includes the optional cooling unit, no additional connections are required. Any condensate produced during cooling operation will be evaporated inside the cabinet and the humid air be exhausted via the ventilation of the cabinet.

2.3.4 Installation of the Grid for Petri Dishes (optionally)

Scope of delivery:


Figure 5: Grid for petri dishes with installation material

- | | |
|--------------------------|-----------|
| 1. Fastening screws | (12 pcs.) |
| 2. Mounting rails | (4 pcs.) |
| 3. Angle bracket | (4 pcs.) |
| 4. Grid for petri dishes | (1 pcs.) |

2.3.4.1 Installation of the Mounting Rails (4 pieces)

1. Attach the first mounting rail with the upper screw in the provided position of the side wall (slightly tighten the screw).

 Care for the correct mounting position (see figure 6)

Figure 6: Mounting position of the mounting rail

2. Insert the remaining 2 screws into the screw holes and tighten them slightly.
3. Securely tighten all screws.
4. Install the remaining 3 mounting rails in the same order.

2.3.4.2 Installation of the Angle Bracket

1. First attach the upper hook of the angle bracket in the recess of the mounting rail at the desired height (see figure 7).


Figure 7: Installation of the angle bracket

2. Now press the angle bracket in such a way that the lower end of it snaps in the mounting rail – one recess notch under it (see figure 8).


Figure 8: Assembled angle bracket

3. Install the remaining 3 angle brackets in the same order at the other mounting rails.

 Take care for the equal height of the angle brackets.

4. Insert the grid for petri dishes into the device.

 The air flow of the device ventilator inside must not be blocked!

 Always mount the tray insert below the device ventilator located on the back of the device.

 Make sure that the grid for petri dishes has been securely mounted in horizontal position.

5. If necessary:
 - remove the grid for petri dishes
 - remove the incorrectly positioned angle bracket by pressing it slightly
 - install the angle bracket in a new position again

2.4 Equipment, Mounting and Loading of a Tray

Sartorius Stedim Biotech GmbH offers trays for the CERTOMAT® Tplus fixed with sets of universal clamps for Erlenmeyer flasks of different sizes, or trays without assemblies. On delivery or after the first set-up of the CERTOMAT® Tplus the tray can be premounted ready for operation. If this includes all clamps and mountings required for first use, you can insert the containers and start operation, as shown in chapter 3.2 Equipment, Mounting and Loading of a Tray. If you need to change the tray or its equipment for your application or if you want to insert another tray heed the following information.

2.4.1 Equipping the Tray

- We recommend to mount all clamps and mounting sets for your containers on the tray before you place the tray in the Tplus. So they can be handled easier. However, the flasks, bottles, etc. with the media (especially heavy containers) should be placed on the tray after it has been attached to the tray support.
- Take care of the dimensions of the incubation chamber, about 890 x 495 x 650 mm (W x H x D)

Mind the weight of the fully equipped tray. It may be very heavy. The shaking containers may hinder the handling. You might jolt with the tray and damage the containers placed on the tray. Media may be released unintentionally and may cause personnel hazards.

1. For a universal tray without clamps, etc., place and screw tight the holders and mountings required. In order to change an existing tray setup, first release all parts which are not needed.
- Clamps for Erlenmeyer flasks are available in different sizes which can be individually combined. For test tubes, bottles, beakers or separation funnels, mounting systems, such as test tube racks or universal mounting sets, are available. Detailed information about available container holders and mounting systems you will find in section 12.2.3
2. When assembling clamps and other devices on the tray, make sure to distribute them evenly. This will avoid uneven loads and ensure that the containers can be safely fixed on the tray / in the clamps.

2.4.2 Adjustment of the Counter Weight to the Intended Load

 Upon delivery the CERTOMAT® Tplus is preadjusted for a load of 10 kg. If the intended load is different from this, the counter weight has to be readjusted correspondingly. Thus the vibrations, which may occur during operation, can be minimized.

1. Weigh the completely equipped tray with all shaking containers attached.
2. Switch off the CERTOMAT® Tplus at the main switch.
3. To adjust the counter weight with the variable tool for the load open the front flap under the glass doors downwards. The front flap is fixed with magnets. The variable tool is located behind this flap.
4. Turn the shaking table by hand to its rearmost position.
5. Block the drive of the shaking table by putting the delivered hexagon socket screw key through the boring of the shaking table (Fig. 9 and 10, pos. 4). Carefully turn the shaking table until the drive has totally blocked.


Fig. 9:: Adjustment of the counter weight corresponding to the load

6. Put the variable tool (1) through the opening under the front flap in the housing and guide it on the axis of the counter weight (2);
7. Turn the counter weight with the variable tool in the clockwise direction (for load <10 kg) or turn it out anticlockwise (load >10kg).
 - The adjustment of the counter weight which is necessary for a certain load is marked on the variable tool.
8. Remove the hexagon socket screw key and the variable tool and close the front flap.
9. Mount the loaded tray.
10. Switch on the CERTOMAT® Tplus.

 The CERTOMAT® Tplus starts only if both doors are closed.

11. Test for the silent operation at your intended shaking speed.
12. If a re-adjustment is necessary, stop the shaker, switch off the CERTOMAT® Tplus, remove the tray and repeat steps 3 – 11.

1. Variable tool
2. Counter weight
3. Axis
4. Hexagon socket screw key

2.4.3 Mounting and Dismounting a Tray


Fig. 10: Mounting of the tray on the shaking table (tray shown without mounting systems and containers)

1. Place the tray onto the shaking table. Push it with the indents (1a/1b) against the rear knurled screws (1c/1d) and the clamp (2a) of the shaking table.

Please mind the weight of the tray if you have already loaded it with shaking containers.

2. Push down (2b) the front edge of the tray so that it properly fits into the front clamp (2c).
3. The borings (3a/3b) will precisely fit onto the locating pins of the tray support (3c/3d) and locate the tray.

Please take care that the locating pins (3c/3d) are placed evenly in the tray. Only then the tray is situated correctly on the shaking table.

If necessary, you can adjust the locating pins as follows:

1. Release the screws of the locating pins (3c/3d) with the hexagon socket screw key.
2. If necessary, move the tray slightly. The locating pins (3c/3d) position themselves in the centre of the borings (3a/3b) of the tray.
3. Tighten the screws of the locating pins (3c/3d) carefully again after adjusting the tolerance.
4. In order to improve the safe position of the tray – especially at high loads in connection with high speed – the provided flat headed screws M5 x 8 must be screwed in the prepared borings.

2.5 Mounting of Trays and Mounting Sets

2.5.1 Assembly and Loading of Finally Equipped Sets of Trays


Fig 11: Tray type F with clamps for Erlenmeyer flasks

The trays type E and F are available as completely equipped trays. Detailed information you will find in the "Ordering Information" in the supplement.

1. Place the tray in the CERTOMAT® Tplus and fasten it as explained in chapter 2.4.3 Mounting and Dismounting a Tray.
2. Insert the Erlenmeyer flasks in the clamps (you can also use bottles, beakers etc. which fit in the clamps).

Please take care of the distribution of the loads and distribute them evenly starting from the centre of the tray outward.

Before starting the shaker check the proper fixation of all containers. They must not get loose and flung off during operation of the shaker.

2.5.2 Mounting of Test Tube Racks


Fig 12: Test tube rack assembly

1. Loosen the clamping screw of the test tube rack and turn the test tube holder, until the screw-holes for attachment of the rack to the tray are accessible. If necessary, you may completely remove the holder from the rack.
2. Place the test tube rack onto the universal tray and tighten with the mounting screws. Consider the number of racks required for your application and distribute the racks evenly on the tray starting from the centre of the tray outward.
3. As far as being disassembled, mount again the test tube holders to their racks. Turn the holders into their working position (usually vertical). Locate with their clamping screw.
4. The test tube holders can be fixed at any angle of use as required. Make sure that no liquid escapes from the vessels during operation.

Avoid contact of the vessels with one another (danger of breakage!).

2.5.3 Mounting of the Universal Mounting System

The universal mounting set is available for the trays of size EU and FU. Enclosed you will find further information about the component parts. Basic parts are two side frames which are adapted to the corresponding tray and the rods for connection of the side frames:

- Set of side parts type B-2, art. no. 8854238, for tray EU
- Set of side parts type B-3, art. no. 8854243, for tray FU

Universal clamping rods and mounting sets for separation funnels serve as mounting systems

- Universal clamping rods type U, art.-no. 8854254, for side part set B-2 / B-3
- Mounting set type S-1 for separating funnels 50 and 100 ml, art. no. 8854262, with clamps and retaining springs and 1 universal clamping rod type U
- Mounting set type S-1 for separating funnels 250, 500 and 1000 ml, art. no. 8854270, with clamps and retaining springs and 1 universal clamping rod type U

1. Place the side parts (1) onto the tray and screw them onto the tray with the mounting screws as shown at pos. (1a). For each side part 4 screws are supplied to provide proper tightening. Then mount connecting rods as shown at pos. (1b).
2. For a mounting system with universal clamping rods attach the clamping rods at the side frame rods as shown at pos. 3. You'll need at least two rods and you can add as many as required for placement of your vessels.
3. Place the shaking containers between the clamping rods. Distribute the containers evenly on the tray, starting from the centre outward. You can loosen the clamping screws 3a/3b of the rods to readjust the rod and properly position the flasks, vessels and beakers, etc. Then carefully tighten the clamping screws (3a/3b) so that the containers do not loosen during operation.
4. For the mounting set for separation funnels attach one universal clamping rod as shown above. It will serve as support. Then mount the retaining springs (2a) to the clamping rod for separation funnels (2). Depending on the size of the separation funnels you can mount up to 5 retaining springs to the clamping rod (pos. 2a).
5. Place the separation funnels onto the mounting set. You can readjust the clamping rods as required using their clamping screws attached.


Fig. 13: Assembly of universal rods and mounting sets for separating funnels

3 Operating Information

3.1 Safety Notes

Media used with incubation shakers like the CERTOMAT® Tplus can be hazardous. This is especially true for caustic, inflammable or infectious media, for instance. Since the risks and the correlated safety requirements depend on the media, they cannot be discussed herein. Your company must issue applicable safety guidelines, distribute them to the laboratory staff and ensure that they are strictly followed.

The following safety instructions must be observed:

When hazardous agents are to be used only authorized personnel should have access to the laboratory or to the working area, respectively, and be allowed to operate the CERTOMAT® Tplus. It must be guaranteed that unauthorized persons do not have access to the unit.

Operators should wear adequate working clothes and personal protection, such as gloves, protective goggles, or breath protection, for instance.

Only equipment and accessories must be used which Sartorius Stedim Biotech GmbH has released for the Tplus. Before starting operation all parts should be checked for damages, especially the glass parts. You must not use any damaged or defective equipment.

Do not start the shaking drive without a tray installed. Always switch-off the shaking drive before mounting or removing the tray or any of the containers to be incubated.

Properly fasten the tray and all mountings and holders to prevent that the containers can loosen during operation. When loading the tray, evenly distribute the containers on the tray (i.e. balanced arrangement) and ensure that they cannot hit against each other.

Reduce the shaking speed as far as necessary, if resonance effects occur and vibrations are transmitted to the working place or neighbouring equipment or if the stack arrangement becomes unstable.

Critical conditions are possible even during normal operation. At specific shaking speeds and depending on the load the action of the shaking drive can cause vibrations due to „self-resonance“ effects. Such vibrations cannot be fully compensated by the internal balance weights of the tray support and will be transmitted to the Tplus and the working place.

If media are incubated in open containers you should lower the shaking speed as far as possible, to avoid unintended spillage of the media. If the media are released unintentionally, i.e. due to breaking of glass vessels, for instance, you should clean the spill tray and incubation chamber. Remove broken glass and any foreign materials as soon as possible and prevent those parts from entering the gap between the tray support and the spill tray or from being introduced into the drive unit. Further information is given in chapter 11.3 Maintenance and Trouble-Shooting.

3.2 Equipment, Mounting and Loading of a Tray

3.2.1 Equipping the Tray

1. If you have changed the mounting place (e. g. cooling room <-> laboratory) please wait until the device has warmed up to room temperature before you start it. Otherwise the air humidity may condensate in the device leading to functional disturbances.

 Following the transport of a CERTOMAT® Tplus with integrated cooling unit, the device may only be put into operation after a waiting period of an hour. This avoids danger to the cooling unit.

2. Connect the CERTOMAT® Tplus with the mains cable to a suited socket.
3. By turning the mains switch (2) the electronics are supplied with line voltage.

 After switching on the control unit the integrated DCU system initialises. This takes approx. 1-2 minutes. The system is ready for operation when the "Main" menu appears on the screen. The "Main" menu gives a schematic overview of the system with symbols for temperature, shaker speed, status of the illumination and cooling, for example, as well as the measured value displays associated with them.

 The "Main" menu is the starting point for operation of the CERTOMAT® Tplus.

Fig. 14: View of the "touch display" and mains switch

1. Touch panel
2. Mains switch

3.2.2 System behaviour at the Start

After switching on the control unit and program start as well as following the return of the power after a power failure the measuring and regulating system is in a defined basic state:

- The system configuration is loaded as it was preset ex works or implemented later.
- User-defined parameters of a previous process can be stored and called up again:
 - Set values
 - Calibration parameters
 - Programming (batches)
- All controllers are switched off ("off").
- Timer sequences (e.g. the order of events of a stored program) are stopped ("off").

In the case of downtimes, the switch-on behaviour depends on the outputs and system functions that have a direct effect on the connected terminal (controller, timer, etc.) and on the type and duration of the downtime. The following types of downtime are differentiated:

- Switching on and off at the "Mains" main switch on the control unit.
- Downtimes caused by failure of the power supply from the laboratory connection (power failure).

A maximum duration or "Failtime" for power failures can be set in the "System Parameters" submenu of the main "Settings" function.


Fig. 15:: "System Parameters" submenu
 → For a description, see chap. 8 "Settings"

When the voltage returns after a power failure that was shorter than "Failtime", the system continues to work as it was:

- The system shows the "Power Failure" error message with the time of the failure.

- Controllers continue to work with the selected set value.
- Timers and set value profiles continue to be worked through.

If the mains power failure persists longer than the selected "Failtime", the DCU system behaves as if the user had switched the device off normally, i.e. it takes the defined basic state.

After the next restart the alarm message "Power Failure, Fermentation stopped" then appears (→ alarm messages in the Appendix) with details of the date and time at which the power failure occurred.

4 Basics of the Operation of the "Touch Display"

4.1 Layout of the main screen

The "Main" screen provides a graphic overview of the installation of the controlled device where symbols indicate the typical arrangement in the shaker.


Fig. 16: "Main" menu, CERTOMAT® Tplus

*) Schematic configuration, functional elements, tags and parameters depend on the configuration.

4.1.1 Header

The header of the screen shows status information only:

Selected main function: – Example "Main"	Main 01: Idle 0:00 h:m 2008-04-02 10:11:08
	Display of the active main function in the working area
Main	Display of the selected main menu
01: Idle	Program display (batches) Status of the program selected via the main "Batches" function
0:00 h:m	Access to the submenu for entering a time (reference time for process start)
2008-04-02 11:28:11	Date in the format [yyyy-mm-dd]; Time in the format [hh:mm:ss]

Fig. 17: Example: Header - "main" screen

4.1.2 Working Area

The "working area" shows the functional elements and submenus of the active main function:

- Preselected process values with their current measured value or set value
- Controllers, e.g. for shaking unit, etc., with current set values
- The current actual values and actual states are displayed in the "main" screen.


Fig. 18: Example: Working area - "Main" screen

- The DCU system is operated directly on the display by selecting a main function and its associated submenus. The functional elements in the working range and the main function keys in the footer contain "touch keys". Touching them activates associated submenus, e.g. for the input of data and set values or the selection of operating modes.
- Available functions, "tag" names, parameters and submenus depend on the configuration of the DCU system.

4.1.3 Footer

The footer contains the "main function keys" for switching between the "main functions":

Selection of "Main" via the main function key


Fig. 19: Example: Footer - "Main" screen

Display form:

- Selected main function: Key is light grey, pressed down
- Non selected functions: Keys are dark grey, raised

- Selection of "Main" via the main function key

4.2 Display of the Functional Elements

The display of the functional elements in the working area indicates their current status and intended use:

Symbol	Display	Meaning, Use
	Key with grey underline	Measured value recording or output of the functional element is inactive, with measured value or control variable as shown
	Key with green underline	Measured value recording or output of the functional element is active, with measured value or control variable as shown
	Key with yellow underline	Measured value recording or output of the functional element is active when "Manual" operating mode is selected; Automatic control is not possible
	TEMP "off" → auto "on" Line grey → green	Direct access from the "Main" menu to the submenu for selection of the operating mode
	TEMP "off" → manual "on" Line grey → yellow	Direct access from main "Controller" menu to the submenu for selection of the operating mode

Display of functional elements

Examples of functional elements, abbreviations, measured values, operating values and submenus that can be called up through selection of the "touch keys" [→ part of the individual main functions].

4.3 Overview of the Main Function Keys

Key, Symbol	Meaning, Use
"Main" function 	Start screen with graphical overview of the controlled device: <ul style="list-style-type: none"> – Display of the components of the current configuration – Overview of the measured variables and process parameters – Direct access to important menus for operating inputs
Main "Controller" function 	Operating and parameterization menus for controllers, e.g.: <ul style="list-style-type: none"> – Temperature control, TEMP – Speed control of the shaker (SHAKER)
Main "Batches" function 	Menu for dealing with the programs <ul style="list-style-type: none"> – Input of programs (sequences) for the automatic operation of the CERTOMAT® Tplus – Start and administration of programs
Main "Settings" function 	Basic system settings such as: <ul style="list-style-type: none"> – Measuring ranges of process values – Manual operation, e.g. for inputs and outputs, controllers, etc. – External communication (e.g. with printers, external computers) – Selection, changing configurations (password-protected, only by authorized service personnel)
„Remote" operation 	Operation with external computer systems (central computer) <ul style="list-style-type: none"> – Pressing the main function key switches over to "Remote operation"; For notes about the configuration, see [→ Chapter 8 Main "Settings" Function]
„Alarm" display 	Overview table of the alarms that occur: <ul style="list-style-type: none"> – If alarms occur, the symbol changes colour and an acoustic signal sounds. – Display red : contains alarms which haven't been acknowledged – Pressing the main function key opens an overview menu of all alarm messages.

Overview and display of the main function keys

All main functions can be selected at any time while a process is running.
 The title of the selected main function displayed in the working area also appears in the header.

5 "Main" Function

Fast access to all controllers (TEMP, SHAKER) and switches (LIGHT, ILLUM, COOL) is obtained via the "Main" function key in the footer. After pressing the key an overview menu appears first.

Then open the controller operating images (TEMP, SHAKER) and the switches (LIGHT, ILLUM, COOL) by pressing the functional element.


Fig. 20: "Main" menu of the CERTOMAT® Tplus

5.1 Functional Elements of the "Main" Menu

Key, Symbol	Meaning, Use
	<ul style="list-style-type: none"> • Measured temperature value in the CERTOMAT® Tplus – Input of the temperature set value – Operating mode selection for TEMP controller – Continue with the controller menu
	<ul style="list-style-type: none"> • Shaker speed [rpm] – Input of the set value for the shaker speed – Operating mode selection for SHAKER controller – Continue with the SHAKER controller menu
	<ul style="list-style-type: none"> • Working light ON / OFF – Manual operation only
	<ul style="list-style-type: none"> • Photo illumination ON / OFF – Can be switched manually – Can be switched by program control
	<ul style="list-style-type: none"> • Cooling ON / OFF – Can be switched manually – Can be switched by program control

5.2 Starting an Incubation from the "Main" Menu

In order to start an incubation from the "Main" menu, proceed as follows.

Beispiel:

- The following operating values are specified:
 - Temperature: 37 °C
 - Speed: 50 rpm.
 - Working light: OFF
 - Photo illumination: ON
 - Cooling: OFF

5.2.1 Entering the Temperature

- Press the "TEMP" key.
The input menu of the temperature controller appears.


Fig. 21: "Controller TEMP" menu

- Enter the temperature (37 °C) via the screen keypad, set the "Mode" to "auto" and confirm with "ok".

The heating mode starts.
 The "TEMP" key shows the green underline for operating mode "auto" and the entered temperature (37 °C) is displayed

5.2.2 Entering the Speed

- Press the "SHAKER" key.
The input menu of the speed controller appears.

Fig. 22: "Controller SHAKER" menu

- Enter speed (50 rpm) via the screen keypad, set the "Mode" to "auto" and confirm with "ok"

SHAKER /min
50

The shaking operation starts. The "SHAKER" key shows the green underline for operating mode "auto" and the entered speed (50 rpm) is displayed.

5.2.3 Switching the Working Light on or off

- Press the "LIGHT" key.
The input menu for the working light appears.

Fig. 23: Entry LIGHT "On"/ "Off"

- Switch off working light by pressing "Off".

LIGHT

The working light is switched off. The "LIGHT" key indicates the grey underline for operating mode "Off".

5.2.4 Switching the Photo Illumination on or off

- Press the "ILLUM" key.
The input menu for the photo illumination (ILLUM) appears.


Fig. 24: Entry for ILLUM "On"/ "Off"

- Switch on photo illumination by pressing "On"

ILLUM

The photo illumination is switched on. The "ILLUM" key shows the green underline for operating mode "On".

5.2.5 Switching the Cooling on or off

- Press the "COOL" key.
The input menu for the cooling appears.



Fig. 25: COOL input "On"/ "Off"

- Switch off cooling by pressing "Off".

COOL

The cooling is switched off. The "COOL" key shows the grey underline for operating mode "Off".

View of the "Main" Menu during the Running Process.

5.3 Editing individual Values of an Incubation from the "Main" Menu

While the incubation is running, changes can be made to individual values (temperature, speed) or states (working light, photo illumination, cooling) at any time.

A running incubation process does not have to be stopped beforehand in order to do so.

Entering "0" rpm does **not** completely stop the shaker! Motor will run at approx. 5% of max. speed.

The changing of values and states is carried out as described in chapter 5.2.1 to 5.2.5

5.4 Ending an incubation from the "Main" menu

In order to end an incubation from the "Main" menu, proceed as follows.

5.4.1 Switching off the heating

- Press the "TEMP" key.
 - The input menu of the temperature controller appears.
- Set the "Mode" to "off" and confirm with "ok".

TEMP °C The heating is switched off.
37.0 The "TEMP" key shows the grey underline for operating mode "off". The entered temperature (37 °C) remains stored.

5.4.2 Switching off the Shaker (Speed 0 rpm)

- Press the "SHAKER" key.
 - The input menu of the speed controller appears.
- Set the "Mode" to "off" and confirm with "ok".

SHAKER /min The shaker drive is switched off.
0 The "SHAKER" key shows the grey underline for operating mode "off".

5.4.3 Switching the Working Light, Photo Illumination and Cooling on or off

The switching states of the working light, photo illumination and cooling are changed as described in chapter 5.2.3 to 5.2.5

6 Main "Controller" Function

Access to all controllers is gained via the main "Controller" function key in the footer.
 After pressing the key an overview menu appears first that shows the available controllers.
 Then open the controller operating images by pressing the functional element of the individual controller.


Fig. 26: Main "Controller" menu of the CERTOMAT® Tplus

6.1 Functional Elements of the Main "Controller" Menu

Key, Symbol	Meaning, Use
	<ul style="list-style-type: none"> • Set value controller for temperature in the CERTOMAT® Tplus – Works with the temperature of the CERTOMAT® Tplus as control input, measured with the sensor that is mounted in the air stream
	<ul style="list-style-type: none"> • Set value controller for CERTOMAT® Tplus – Controller for shaking speed of CERTOMAT® Tplus

A system configuration takes into account suitable operating values and necessary limitations of the controller for the controlled terminal [→ Configuration documentation]. At the same time each controller is parameterized in such a way that it enables stable control.

6.2 Access to the Controller

Example: Temperature Controller

6.2.1 Entering the Set Temperature Value

- Press the "TEMP" key.
 The input menu of the temperature controller appears.


Fig. 27: "Controller TEMP" menu

- Press the "Setpoint" key.
 The "TEMP Setpoint" input menu appears.


Fig. 28: "TEMP Setpoint" menu

- Enter new set value via the screen keypad and confirm with "ok".

6.2.2 Setting the Alarm Limits

- Press the "Alarm Param" key.
The "Alarm TEMP" input menu appears.


Fig. 29: "Alarm TEMP" menu

Pos.	Meaning
1	After pressing the "Delay" key the menu appears for entering the delay time for the "Highlimit" alarm triggering temperature in seconds.
2	Current actual value of the temperature
3	After pressing the "Delay" key the menu appears for entering the delay time for the "Lowlimit" alarm triggering temperature in seconds.
4	Optical display of the alarm status <ul style="list-style-type: none"> Bell appears crossed out <ul style="list-style-type: none"> – Alarm message for monitoring the alarm limits is deactivated Bell appears not crossed out <ul style="list-style-type: none"> – Alarm message for monitoring of the alarm limits is activated
5	Switch for activation (enabled) or deactivation (disabled) of the alarm limit monitoring
6	After pressing the "Lowlimit" key the menu for entering the lower temperature limit in °C appears.
7	After pressing the "Highlimit" key the menu for the input of the upper temperature limit in °C appears.

- After entering the alarm limits activate the alarm by pressing the "Alarm" key (Fig. 29, Pos. 5). Display shows "Alarm enabled". In addition the alarm limits are indicated by red triangles (Fig. 30, Pos. 1 and 2) in the bar graph.


Fig. 30: "Alarm TEMP" menu

Pos.	Meaning
1	Optical identifier (red triangle) of the upper temperature limit
2	Optical identifier (red triangle) of the lower temperature limit
3	Optical display of the alarm status <ul style="list-style-type: none"> – Monitoring of the alarm limits activated.

- Leave the "Alarm TEMP" menu by pressing the "red cross". This returns you to the "Controller Temp" menu (Fig. 27).

6.2.3 Select "Controller Operating Mode".

- Press the "off" key.
The "Controller Mode" input menu appears.


Fig. 31: "Controller Mode" input

- Select the desired operating mode for the controller by pressing the "auto" or "manual" key.
- Leave the subsequent controller menu by pressing the "ok" key.
- In order to display all process parameters return to the "Main" menu with the "Main" key in the footer.

7 Main "Batches" Function


Fig. 32: Main "Batches" menu of the CERTOMAT® Tplus

7.1 Functional Elements of the Main "Batches" Menu

Key, Symbol	Meaning, Use
	<ul style="list-style-type: none"> • Measured temperature value in the CERTOMAT® Tplus – Input of the temperature set value – Operating mode selection for TEMP controller
	<ul style="list-style-type: none"> • Shaker speed [rpm] – Input of the set value for the shaker speed – Operating mode selection for SHAKER controller
	<ul style="list-style-type: none"> • Working light ON / OFF – Manual operation only
	<ul style="list-style-type: none"> • Photo illumination ON / OFF – Can be switched manually – Can be switched by program control
	<ul style="list-style-type: none"> • Cooling ON / OFF – Can be switched manually – Can be switched by program control

7.2 Example of Creating a Program

- A total of 20 different programs can be stored in the device. Each program consists in turn of a "Pre-step" and four further steps that are worked through one after the other. For each step the values for cooling, photo illumination (ILLUM), speed and temperature can be set independently of one another. Each step has a defined running time.
- If only one or two steps are required, the remaining steps can be defined with a running time of 00:00 so that the program stops directly afterwards. A repeat function for endless loops is provided.
- The "Pre-step" can only be switched on or off (no running time). It defines the conditions (cooling, photo illumination (ILLUM), speed and temperature) from pressing the [Start] key until the program begins (Step 1).
- If "Pre-step" is switched off, cooling, photo illumination (ILLUM), speed and temperature are switched off until the program begins.
- The main application area for the "Pre-step" is in the scheduled timely provision of cultures, e.g. after a weekend or the following morning.

Advantages:

- A program only has to be defined once with running times and it can be started any time later. The running times of the program always remain the same, only the duration of the "Pre-step" varies.
- Comfortable starting of night or weekend programs through the day or before the weekend.
- Conditions are reproducible before the start of the program because these are held by the device according to the settings in the "Pre-step".
- Automatic program start without the presence of personnel.

Example:

On Monday morning at 9:00 o'clock any cultures have to be ready. The program (Step 1 ... Step 4) for the production of these cultures needs altogether 36 h. So the program must be started manually without use of the "Pre-step" on Saturday evening at 7:00 o'clock.

The "Pre-step" makes it possible to choose any time for the start and to keep the material cooled under the desired conditions.

Programm 01

	Speed/ Duration [h]	Temperature [°C]	Drehzahl [1/min]	Light (ILLUM)	Cooling
Pre-Step	on	20	0	off	on
Step 1	7	30	200	on	on
Step 2	10	40	100	off	off
Step 3	20	10	60	on	on
Step 4	1	10	200	on	on








Fig. 33: Example for programming


7.3 Entering Program Data

Step	Key, Symbol, Display	Input
Call up main "Batches" menu		Press the "Batches" key.
		
Select program		Press the "Program" key.
		
Enter program no. via screen keypad and confirm.		Enter program no. via screen keypad and confirm with "ok" key.
		

Step	Key, Symbol, Display	Input
Select mode		Press the "Mode" key.


Enter program data		Press the "Modify" key.
--------------------	---	-------------------------

Define pre-step










Enter speed during the pre-step		Press the "Value" key in "Pre-Step" line / "SHAKER" column.
---------------------------------	---	---


Enter speed via screen keypad.		Enter speed via screen numeric keypad and confirm with "ok" key.
--------------------------------	---	--

If necessary, switch alarm message on or off	 → deactivated  → activated	Press "Alarm" key for desired status Note: The "Highlimit" and "Lowlimit" alarm limits are specified by the corresponding settings of the controller. These can be set in the "Main" menu or main "Controller" menu
--	--	--

 The alarm limits "Highlimit" and "Lowlimit" are specified by the corresponding settings of the controller. These can be set in the "Main" menu or main "Controller" menu.

Step	Key, Symbol, Display	Input
Enter temperature during the pre-step		Press "Value" key in "Pre-Step" line / "TEMP" column.
		
Enter temperature via numeric keypad on the screen.		Enter temperature via screen keypad and confirm with "ok" key.
Switch alarm message on or off	 → deactivated  → activated	Press "Alarm" key for desired status.
<p> The "High limit" and "Low limit" alarm limits are specified by the corresponding settings of the controller. These can be set in the "Main" menu or main "Controller" menu.</p>		
Switch cooling on or off during the pre-step		Press the "Value" key in the "Pre-Step" line / "COOL" column.
		
Switch cooling on or off.	 	Press "On" or "Off" key for desired status.
Switch photo illumination on or off during the pre-step		Press "Value" key in "Pre-Step" line / "ILLUM" column.
		
Switch photo illumination on or off.	 	Press "On" or "Off" key for desired status.

Step	Key, Symbol, Display	Input				
Define step 1						
Enter operating time of step 1	<table border="1"> <tr><td>Duration</td></tr> <tr><td>0 Pre-Step</td></tr> <tr><td>1 Value</td></tr> <tr><td>0d 00:30</td></tr> </table>	Duration	0 Pre-Step	1 Value	0d 00:30	Press "Value" key in "Step 1" line / "Duration" column.
Duration						
0 Pre-Step						
1 Value						
0d 00:30						

Enter operating time [hh:mm] via screen keypad.		Enter operating time [hh:mm] via screen keypad and confirm with "ok" key.
---	--	---

Input of the duration is basically made in hours and minutes [hh:mm]. After input, operating times greater than 24 h are automatically converted into days, hours and minutes.


The remaining values (speed, temperature, cooling and photo illumination) for the programming of step 1 must be entered in the same way as the programming of the pre-step.

Define steps 2, 3 and 4

The values for the programming of steps 2, 3 and 4 must be entered in the same way as the programming of the pre-step or step 1.







If steps 2, 3 and 4 are not required, the operating time of the respective step must be set to "0".


Leave "Batch Parameters" menu.		Press "Close" key.
--------------------------------	--	--------------------

7.4 Specifying the Start Time of a Program


7.4.1 Starting the Program immediately

Step	Key, Symbol, Display	Input
Call up main "Batches" menu.		Press "Batches" key.
		
Select mode.		Press "Mode" key.
		
Start the program.		Press "Start" key.
		
Enter number of cycles.		Press "Cycles" key.
		
Enter number of cycles via screen keypad.	 → 	Enter number of cycles via screen keypad and confirm with "ok" key.


Start the program.  Press the "Start" key.

- The process starts immediately.


7.4.2 Starting the Program at a specified Time


Step	Key, Symbol, Display	Input
Call up main "Batches" menu.		Press the "Batches" key.


Select mode.  Press the "Mode" key.


Call up the "Start at time" menu.


Press the "Start at time" key.


Enter start time.


Press the "Time" key.


Enter start time in hours, minutes, seconds [hh:mm:ss] via screen keypad.

Enter start time in hours, minutes, seconds [hh:mm:ss] via screen keypad and confirm with "ok" key.

Enter start date.

Press the "Date" key.


Step

Key, Symbol, Display


Input


Enter start date as year, month, day [yyyy-mm-ss] via screen keypad.


Enter start date as year, month, day [yyyy-mm-ss] via screen keypad and confirm with "ok" key.

Enter number of cycles

Press "Cycles" key.

Enter number of cycles via screen keypad.


Enter number of cycles via screen keypad and confirm with "ok" key.

Start the program.





Press the "Start" key.

- The process starts at the specified time. The conditions from the "pre-step" are carried out until then.

7.5 Stopping / Cancelling a Program

To stop or cancel a running program, proceed as follows:







Step	Key, Symbol, Display	Input
Call up main "Batches" menu.		Press the "Batches" key.
		
Select mode.		Press the "Mode" key.
		
Stop program.		Press the "Stop" key.
		


-  Only a program stop/cancellation takes place. All process values/states active at the time of the program stop/cancellation (temperature, speed, photo illumination and cooling) remain active. The shaker continues to run with the values entered before the program stop/cancellation.
-  Changes are carried out directly in the respective controller for temperature (TEMP) or speed (SHAKER) or at the photo illumination (ILLUM) or cooling (COOL) switch.

- Carry out change of temperature or speed by pressing the "TEMP" or "SHAKER" key.
 - Enter new value and confirm with "ok" or set controller mode to "off" and with confirm "ok".
- Make changes to the photo illumination or cooling by pressing the "ILLUM" or "COOL" key.
 - Switch photo illumination or cooling off or on by pressing the "Off" or "On" key.

7.6 Editing the Program

In order to edit a stored program, proceed as follows:

Step	Key, Symbol, Display	Input
Call up main "Batches" menu.		Press the "Batches" key.
		
Select program to be edited.		Press the "Program" key.
		
Enter program no. via screen keypad and confirm.		Enter program no. via screen keypad and confirm with "ok" key.
		
Select mode.		Press the "Mode" key.
		
Edit program data		Press the "Modify" key.

- Adjust values according to requirements.
 - Edit values as described in section 7.3 "Entering Program Data".

Leave "Batch Parameters" menu.

Press "Close" key.

8 Main "Settings" Function

-  The main funktion "Settings" enables interventions in the system configuration. Malfunctions with unforeseeable effects on safe operation can arise from settings that are not permissible or unsuitable for a specific end terminal.
-  Settings that have an effect on safe operation are protected with passwords. They are only allowed to be changed by experienced or specially trained personnel.
-  The standard password [→ Appendix] may only be passed on to authorized users and the service password [→ separate notification] may only be passed on to authorized service employees and administrators.

8.1 General Information

The DCU system provides various functions for system maintenance and the elimination of faults via the main function system settings ("Settings"):

- General settings such as date, time, fault waiting time or "Failtime", parameterization of the communication with external devices ("Internet Configuration").
- Specification of process values ("PV") and their ranges or limits.
- Manual operation, e.g. of digital and analog inputs and outputs or controllers, for the simulation.
- Service function, e.g. for system reset or for selection of the system configuration, insofar as multiple configurations are implemented or enabled.





Fig. 34: Main "Settings" menu of the CERTOMAT® Tplus

8.2 Functional Elements of the Main "Settings" Menu

Key, Symbol	Meaning, Use
	<p>The system settings ("System Parameters") enable general settings on the DCU system.</p> <p> System parameters: Editing of the time and date values is only possible up to 10 min. after system start</p>
	<p>Setting of process value ranges</p> <p> Only authorized personnel (service technicians) are allowed to make settings in this menu.</p>
	<p>Calibration of the temperature sensor</p> <p> Only authorized personnel (service technicians) are allowed to make settings in this menu.</p>
	<p>Access to the configuration menu of the DCU system supplied (e.g. "Factory Reset", "Boot Configuration", "Reboot")</p> <p> Only authorized personnel (service technicians) are allowed to make settings in this menu.</p>
	<p>For function tests during commissioning and for troubleshooting faults, the digital and analog process inputs and outputs, controller and digital counters can be switched to manual operating mode ("Manual Operation").</p> <p> Settings in manual mode have the highest priority; they have priority over the other functions on the associated control elements of the DCU system. Changes to settings are to be carried out solely test purposes upon consultation with the technical service department.</p>

8.3 Operating Menus and Input Window

8.3.1 "System Parameters"

After pressing the "System Parameters" key and subsequent entry of the standard password the following settings can be carried out:

Field	Value	Function
Time	hh:mm:ss	Entry of current time, in "hours:minutes:seconds" format
Date	yyyy-mm-dd	Entry of current date, in "year-month-day" format
Failtime	hh:mm:ss	Entry of power failure time for the system behaviour when it is switched on again
Internet Config	IP-Adresse	Configuration of the Ethernet interface for the connection of external DP facilities, e.g. via the LAN into the laboratory
Beeper	enable/ disable	Switch-on or switch-off of the acoustic alarm signal, affects all alarms

8.3.2 "PV Ranges"

 Only authorized personnel (service technicians) are allowed to carry out settings in this menu.

8.3.3 "Calib. Temperature"

 Only authorized personnel (service technicians) are allowed to carry out settings in this menu.

8.3.4 "Service"

 Only authorized personnel (service technicians) are allowed to carry out settings in this menu.

8.3.5 „Manual Operation“

 Changes to settings may only be carried out for test purposes in agreement with the technical service department. Settings in manual mode have the highest priority, they have priority over the other functions of the associated control elements of the DCU system.

9 "Remote Function"

The remote mode is activated or deactivated by pressing the "Remote" key.


Fig. 35: "Remote" function of the CERTOMAT® Tplus

9.1 Functional Elements of the "Remote" Function

Key, Symbol	Meaning, Use
	Remote mode is deactivated. There is no communication between the external control unit and DCU4 control unit. The remote control is deactivated. The operator cannot access the DCU4 control unit from an external control unit (e.g. MFCS/win).
	Remote mode is activated. Communication between the external control unit and DCU4 control unit takes place. The bidirectional control unit is activated. Remote control of the DCU4 control unit is enabled.

10 Main "Alarms" Function

If an alarm occurs, a message window appears on the current display with details of the date, time and type of error that occurred. The colour of the alarm bell in the footer changes from "white" to "red".


Fig. 36: Alarm message if an alarm occurs

10.1 Functional Elements of the Main "Alarms" Menu

Key, Symbol	Meaning, Use
	No alarm occurred. Any alarms that occurred have been confirmed and/or deleted.
	At least one alarm occurred. Alarm has not yet been confirmed.

By pressing the "Alarm" key the display in the main window switches to the display of alarms. All alarms that occurred during the running process are displayed here


Fig. 37: List of alarms that occurred with a short description of each alarm

10.2 Confirming and/or deleting Alarms

The operator has the following options for reacting to alarm messages:

- Confirm the alarm in the superimposed window with "Acknowledge". The alarm message in the header and red highlighting of the alarm symbol are deleted. The alarm is listed as a "confirmed" alarm ("ACK") in the list of alarms.
- Close the superimposed window without confirming with . The alarm message remains in the header and red highlighting of the alarm symbol is retained. The alarm is listed as an "unconfirmed" alarm ("UNACK") in the list of alarms.
- Call up the alarm list with the main function key . Unconfirmed alarms can be confirmed there with "ACK". Pressing "RST" ("Reset") deletes the message from the alarm list.
- How to react to alarms and other messages depends on their type and cause. The following tables give an overview. If you have any queries about alarms and suitable measures for their remedy, please contact Sartorius Stedim Biotech GmbH.

10.3 Overview of Process Alarms

Text from alarm line	Meaning	Remedy
"Name State alarm ..."	Alarm of digital input with (name)	Confirm alarm with 'ACK' (acknowledge)
„Name low alarm ..."	The process value (name) has undercut its lower alarm limit	Confirm alarm with 'ACK' (acknowledge)
„Name high alarm ..."	The process value (name) has exceed its upper alarm limit	Confirm alarm with 'ACK' (acknowledge)
„Fuse blown ..."	(Automatic) Circuit breaker in the control unit has triggered	Replace defective fuse, switch on automatic circuit breaker or device, inform the service department if it triggers again
„Heater failure ..."	Heating mat or cable defective, overheating protection has triggered	Inform the service department
„TEMP Sens. failure"	Pt-1000 not correctly connected or defective	Check Pt-1000 and connection, replace sensor if necessary

10.4 Process Messages

Text from alarm line	Meaning	Remedy
„Shut down" ... „Name"	"EMERGENCY STOP" key was activated	Switch on device with the designated "name" again (via EMERGENCY STOP key or main switch)
„Name <3mA alarm ..." „(Name of an analog input) < 3 mA alarm"	Analog input with designated "name" has defective cable	Check if cable is trapped on peripheral devices and check that it functions.

10.5 Abbreviations

All designations of the functional elements in the operating menus depend on the configuration of the DCU system. Standard configurations use the following designations:

10.5.1 Process Values

Abbreviation (TAG)	Meaning
TEMP	Temperature on the sensor of the heating
SHAKER	Shaking speed (speed of the shaker for CERTOMAT® systems)

10.5.2 Digital Inputs

Abbreviation (TAG)	Meaning (and measure if necessary)
TEMPC	Warning message, check required of the Pt 1000 in the air stream
HEATC	Alarm signal, short circuit of the heating
SHAKERC	Alarm signal, short circuit of the shaker drive
FUSEC	Alarm signal, safety fuse actuated/replace fuse
DCUON	Status signal from DCU main switch - DCU switched on
SUPON	Status signal from basic device - basic device switched on

10.5.3 Digital Outputs

Abbreviation (TAG)	Meaning (and measure if necessary)
HEAT	Status signal of the heating, "Heating" in operation
TMPON	Heating switched on
SHAON	Shaker switched on

11 Appendix

11.1 Operation of the Internal Cooling Unit

- If the CERTOMAT® Tplus is equipped with the optional cooling unit, you need not care for condensate produced during cooling operation. Existing condensate will be evaporated internally in the cabinet.
- To avoid damages of the refrigerator an integrated overheat protection permanently switches off the unit if the temperature in the incubation room exceeds +38.5 °C.

11.2 Connection of External Equipment

At the rear side of the device the CERTOMAT® Tplus has a 15-pin Sub-D socket for analog signal outputs and a collective alarm contact.

11.2.1 Analog Exits (ANALOG OUT)

- For the external registration of speeds and temperatures, e. g. with a recorder, the CERTOMAT® Tplus has analog signal exits (0 ... 10 V).

The conditions are as follows:

Speed:	0 V corresponds to 0 rpm.
	10 V corresponds to 400 rpm.
Temperature:	0 V corresponds to 0 °C
	10 V corresponds to 100 °C

The moving load at the voltage entries must not be lower than 10 kΩ.

11.2.2 Collective Alarm

Via the 15-pin sub-D "ANALOG OUT" socket two potential-free contacts are available on PIN 7/8/15 which can be used to externally register alarm situations.


Fig. 38: Device rear side CERTOMAT® Tplus

1. Supply entry socket with supply main and shock-proof plug
2. 15-pin sub-D socket for analog signal outputs
 - Pin 1 Signal + / speed
 - Pin 9 Gnd / speed
 - Pin 2 signal + /temperature
 - Pin 10 Gnd / temperature
 - Pin 7 collective alarm NO
 - Pin 8 collective alarm NC
 - Pin 15 collective alarm C
 - Pin 3 signal + /temperature 2
 - Pin 11 Gnd / temperature 2
3. Ethernet connector

11.3 Maintenance and Trouble-Shooting

11.3.1 Electrical Malfunctions

- If the CERTOMAT® Tplus does not start, check the mains supply. The mains cable should be properly connected and the laboratory mains be operative. If other devices in the laboratory are malfunctioning at this time, a mains failure (such as voltage fluctuations) can be the cause.

Check the fuses. They are located in an insert inside the cabinet, behind the right side panel. Disconnect the unit from mains before you open the side doors. A burned-off body or a melted wire will indicate a defective fuse. For replacement note the correct voltage/current of the fuses.

- If the CERTOMAT® Tplus cannot be started although the mains supply is operative and the fuses are in good condition, contact your service representative of Sartorius AG. You will find the telephone number on the next page.

11.3.2 Malfunctions of the Shaking Drive

- If the shaking drive runs noisy or irregularly the belt can be dirty. If the tray support does not turn although the shaking drive is working (you can hear the motor run), the drive belt may be damaged. If the motor is inoperative although the power supply and all necessary adjustments are available, the motor may be damaged.
1. For checking the belt drive, the tray and shaking table must be disassembled. Unscrew the 4 mounting screws in the centre of the shaking table and remove it. You can check the belt through the opening in the cover of the drive.
 2. For the mounting of the shaking table carefully tighten the screws crosswise.
- If the drive belt is damaged (torn) or the motor is suspected to be defective, please contact your service representative of Sartorius AG.

11.3.3 Malfunctions of the Optional Cooling Unit

- The compressor of the cooling unit includes an internal winding protection for overheat detection. This device cuts off the power supply at overheating of the winding. The cooling unit will automatically be enabled for restart after cooling down of the winding. This may last up to 1 hour.
- If the overheat protection cuts-off the cooling unit often or repeatedly, you should contact your service representative of Sartorius AG.

For devices with integrated cooling unit (UHK versions), the following applies for incubations below the ambient temperature:

If the machine is operated continuously for a long period in conditions of relatively high humidity, ice can form on the vaporiser.

The formation of ice on the vaporiser reduces the air circulation.

This can activate the overheating safety cut-out and switch off the machine.

Allow the machine to cool down before turning it on again.

11.3.4 Handling and Eliminating Faults in the DCU System

If technical problems occur, Sartorius Stedim Biotech GmbH requires a qualified report from the user or from the customer service department.

The user or the customer service department can use the "Function bug handling" form (available on request) to describe the problem. As a result, Sartorius Stedim Biotech GmbH can mostly quickly provide the user with a response on how to resolve the problem.

The user should usually register problems to Sartorius Stedim Biotech GmbH via the responsible service department. If contact via the responsible service department is not available, the user can also send the filled out form directly to:

Sartorius Stedim Biotech GmbH
August-Spindler-Straße 11
37079 Göttingen, Germany
Telephone +49.551.308.3118
Fax +49.551.308.3918
info@sartorius-stedim.com
www.sartorius-stedim.com

11.4 Maintenance

Except for necessary cleaning, the replacement of defective fuses and checking of the shaking drive belt the CERTOMAT® Tplus is maintenance-free. If the optional illumination unit is installed, you can replace defective lamps. All bearings of the shaking drive are lubricated for their lifetime. The drive motor is a maintenance-free brushless motor of external rotor type.

Disconnect the mains cable from laboratory grounded socket before cleaning and/or maintenance and service. Unintended starting of the shaking drive during cleaning and maintenance can cause personnel hazards.

- Maintenance and service of the shaking drive, such as the replacement of the drive belt, service of the electrical and electronic equipment or of the optional cooling devices must only be done by qualified and authorized service personnel.

11.4.1 Maintenance of the Illumination Unit

- Maintenance is confined to the checking for and replacement of defective fluorescent lamps and device fuses. If other malfunctions occur, please contact your service representative of Sartorius AG.
1. Switch-off the illumination unit and disconnect from its mains supply.
 2. Pull the defective fluorescent tube out of its holder. Insert the new tube. On delivery commercially available fluorescent tubes "L18W/77 (Fluora)", or fluorescent tubes according to the customer's specifications are installed.
 3. Replace defective fuses, if necessary. The fuses are placed in a fuse assembly behind the door at the right hand side panel. Note to use the same type/size for replacement, as delivered from factory side.

11.5 Cleaning

1. The interior incubation chamber consists of stainless steel. For cleaning, we recommend usual household cleaning agents or alcohol. Do not use aggressive cleaning agents. Chlorous (hypochlorite) agents, for instance, can cause corrosion.
2. Cleaning the spill tray:
 - Switch off the CERTOMAT® Tplus and pull mains plug.
 - Take out the tray.
 - Loosen the 4 screws of the shaking table.
 - Take off the shaking table.
 - The spill tray is now fully accessibly and can be taken out for thorough cleaning.
3. Put the clean spill tray back in place. Mount the shaking table and fix it with 4 screws. Insert the tray. Reconnect the CERTOMAT® Tplus to mains.
4. After incubation of living microorganisms and cells the CERTOMAT® Tplus may need to be disinfected. For disinfection we recommend disinfectants based on alcohols, such as B. Braun Medical AG's Meliseptol®. After cleaning and disinfection with such agents you should thoroughly dry and ventilate the incubation room.

You should remove broken glass, spilled liquids and other impurities immediately after they occur. Please take care that impurities do not dry.

5. Take care that no foreign materials can enter the air-circulation system. Regularly remove dust from the outer covering grids of the built-in ventilator and of the cooling unit.

11.6 Warranty Regulations and Service

Sartorius Stedim Biotech GmbH warrants its products according to the terms specified in the "General Terms and Conditions of Business" and unless other terms were agreed upon in writing. Date of reference is the date of delivery.

- The guarantee includes defects and malfunctioning due to construction, production and material faults, however it does not include defects due to wrong operation and improper treatment and use. Furthermore it does not apply for parts which became damaged due to normal wear and also not for consumer materials.
- The guarantee will lapse if the user or unauthorized third persons modify the device technically or if they use equipment, components and accessories, which have not been passed by Sartorius Stedim Biotech GmbH for the application in connection with the shaking cabinet CERTOMAT® Tplus.

Repairs can be carried out by the Sartorius AG service representative or by authorized personnel of a customer-owned workshop. Please inform your local Sartorius AG sales representative or contact:

Sartorius AG
Servicezentrum Nord
Weender Landstraße 94-108
37075 Göttingen, Germany
Telephone +49.551.308.3729 / 3740 / 3741
Telefax +49.551.308.3730

- Only spare parts passed by Sartorius Stedim Biotech GmbH must be used.
- Defective units can be returned to Sartorius AG.

 Any units to be sent to us for repair have to be in perfect hygienic condition, clean and well packed. If parts were polluted by media and media components, they have to be cleaned, decontaminated, disinfected or sterilised, according to the valid safety rules of this particular application, e.g. for chemical or biological safety.

 The sender has to prove that all safety guidelines have been observed. For this purpose you have to fill in the enclosed decontamination declaration and attach it to the apparatus. No apparatus will be repaired without declaration of decontamination or description of the measures taken.

 Transport damages or subsequent cleaning or disinfection, if required, will be borne by the sender.

11.7 Password System

 Only make this information available to authorized users and to the service department. If necessary, remove this page from the manual and keep it separately.

Certain system functions and settings that should only be accessible to authorized personnel are protected by the standard password system. This includes, for example,

- the setting of the control parameters (e.g. PID) in the control menus,
- in the main "Settings" function
 - the setting of the "PV" process values
 - the setting of the interface parameters for digital and analog process inputs and outputs or the setting of controllers for the simulation in the "Manual Operation" level.

The "Service" submenu of the main "Settings" function can only be accessed with a special service password. This is only available to authorized service personnel.

If a password-protected function is selected, a keypad with the prompt to enter the password appears automatically. The following passwords can be specified:

- Standard password, preset ex works: "[19]".
- Customer-specific standard password: "[____]**).
- Service password: "[_____]" *)

*) These details are sent to you by post or with the [->Technical Documentation]

11.8 Installation Plan, Interface Configuration

An installation plan of the CERTOMAT® Tplus is given in the Appendix of this manual. Additional copies of this operating manual are available on request from:

Sartorius Stedim Biotech GmbH
August-Spindler-Straße 11
37079 Göttingen
Telefon +49.551.308.3118
Fax +49.551.308.3918
info@sartorius-stedim.com
www.sartorius-stedim.com

12 Technical Data, Ordering Information

12.1 Technical Data

12.1.1 Mechanical Design

Shaking cabinet	outside housing 1265 x 710 x 775 mm (W x H x D); incubation room 890 x 495 x 650 mm (W x H x D)
Weight (with equipment)	approx. 230 kg
Tray support	440 x 420 mm (W x D)
Materials	outside cabinet : sheet steel case with acidproof varnish; incubation room inside: stainless steel case
Trays, type/size	Type E/EU (420 x 420 mm) ¹⁾ ; type F/FU : (800 x 420 mm) ¹⁾
Load	depending on tray and mounting system of the Sartorius Stedim Biotech GmbH accessory program

12.1.2 Mains Supply

Laboratory connector	laboratory grounded socket (specific connectors on request)
Mains specifications	230 V (±5 %) 50 Hz or 115 V (±5 %) 60 Hz
Power consumption	about 1.035 kW (including refrigeration)
Photo illumination	approx. 0.39 kW
Working light	approx. 0.0023 kW
Fuses	T6,3A for 230 V, T10A for 115 V
Radio interfer. suppression	according to EN 61326-1:2006

12.1.3 Operating Data

Orbital motion	Ø 25 mm or Ø 50 mm, depending shaker drive	
Shaking speed	40 ... 400 min ⁻¹	
Rotation speed accuracy	max. ±1% (of final value)	
Incubation temperature	heating version	RT +8 °C ... +70 °C
	heating and cooling version	10 °C under RT .. +70 °C
Accuracy of incubation temperature	± 0.1 °C	
Ambient temperature	+10 ... +40 °C (+10 ... + 30 °C at UHK version)	
Ambient humidity	10 ... 60 %; Incubation room 10 ... 80 %, not condensing	

12.1.4 Analog Out

Socket, 15-poles Sub-D	0 ... 10 V, description see chapter 11.2.1
Exactness	± 0,3 V

12.1.5 Collective Alarm

Collective alarm	Contact without potential (directional contact), max. 230 VAC (6.5 A Ohm load)
------------------	--

12.1.6 Option Cooling Unit

Refrigerant	R 134a
Filling	270g

RT = room temperature
(We reserve the right to carry out dimensional and construction modifications)

¹⁾ Trays of modified design for actual locking; rebuilding of old trays for this locking on request.
Please contact the Sartorius AG customer service.

12.2 Ordering Information

12.2.1 Configuration of the Shaking Cabinets

Cat.-no.	Features and Specifications
	Incubation shaker CERTOMAT® Tplus / 25 mm
	Version 230 V / 50 Hz
8865809	Heating version (UH)
8865906	Heating and cooling version (UHK)
	Incubation shaker CERTOMAT® Tplus / 50 mm
	Version 230 V / 50 Hz
8865825	Heating version (UH)
8865922	Heating and cooling version (UHK)

12.2.2 Optional Equipment


Cat.-no.	Features and Specifications
	Support frame
8864489	welded metal sectional-frame construction for up to two CERTOMAT® Tplus shock-resistant stove enamel finish dimensions: 1150 x 220 x 700 mm (W x H x D) 4 height adjustable stands
	Illumination unit for CERTOMAT® Tplus
	consisting of five fluorescent lamps, à 18 W, can be switched off separately, programmable
8861455	Version 230 V / 50 Hz – in connection with cooling UHK
on request	Version 115 V / 60 Hz – can be integrated afterwards
	Darking plates
8861412	Darking plates, set of 2 pieces, excludes ambient light from the incubation chamber, including magnetic strips for mounting
	Grid for petri dishes
8861471	Stainless steel 640 x 880 mm For Petri dishes, height adjustable, with rails, fastening screws and holding device for grid for petri dishes

12.2.3 Accessories

Cat.-no.	Features and Specifications
	Trays with fixed clamps (with clamps for Erlenmeyer flasks)
	Tray E 420 x 420 mm – clamps made of stainless steel
8853533	Tray incl. 39 clamps for Erlenmeyer flasks 100 ml
8853568	Tray incl. 20 clamps for Erlenmeyer flasks 250 ml
8853584	Tray incl. 14 clamps for Erlenmeyer flasks 500 ml
8853606	Tray incl. 9 clamps for Erlenmeyer flasks 1.000 ml
	Tray E 420 x 420 mm – plastic clamps
8853688	Tray incl. 39 plastic clamps for Erlenmeyer flasks 100 ml
8853666	Tray incl. 19 plastic clamps for Erlenmeyer flasks 250 ml
8853677	Tray incl. 14 plastic clamps for Erlenmeyer flasks 500 ml
	Tray F 800 x 420 mm – clamps made of stainless steel
8853738	Tray incl. 74 clamps for Erlenmeyer flasks 100 ml
8853762	Tray incl. 40 clamps for Erlenmeyer flasks 250 ml
8853789	Tray incl. 26 clamps for Erlenmeyer flasks 500 ml
8853800	Tray incl. 15 clamps for Erlenmeyer flasks 1.000 ml
	Universal – Trays
8853002	type EU (420 x 420 mm)
8853037	type FU (800 x 420 mm)
	Sticky tape for universal trays
8860416	Sticky tape "Premium" for universal trays, 30 x 1 mm, roll of 10 m, durable quality
8864497	Sticky tape for universal trays, size 30 x 1 mm, roll of 50 m for fixing of Erlenmeyer flasks, beakers, etc., only for use at low shaking speeds
	Anti-skid layer for universal trays
8864470	Anti-skid layer size 380 x 420 mm for placement of Erlenmeyer flasks, beakers, etc., only for use at low shaking speeds and temperatures up to 50 °C

Test Tube Racks

Fig. 39:
Test tube rack

8853134	Test tube rack for 64 test tubes Ø 14 mm
8853142	Test tube rack for 42 test tubes Ø 16 mm
8853150	Test tube rack for 36 test tubes Ø 18 mm
8853169	Test tube rack for 33 test tubes Ø 20 mm
8853185	Test tube rack for 18 test tubes Ø 25 mm
8853177	Test tube rack for 12 test tubes Ø 30 mm

Racks for centrifuge tubes, with screw cap

8853088	Rack for 42 centrifuge tubes Ø 16 mm
8853096	Rack for 36 centrifuge tubes Ø 18 mm
8853193	Rack for 33 centrifuge tubes Ø 20 mm
8853240	Rack for 12 centrifuge tubes Ø 30 mm

Clamps for Erlenmeyer flasks (for optional mounting onto the universal trays)

8854505	Clamps for Erlenmeyer flasks 25 ml required / possible max. number of clamps for tray type EU - max. 49 pcs. for tray type FU - max. 98 pcs.
8854513	Clamps for Erlenmeyer flasks 50 ml required / possible max. number of clamps for tray type EU - max. 48 pcs. for tray type FU - max. 96 pcs.
8854521	Clamps for Erlenmeyer-flasks 100 ml required / possible max. number of clamps for tray type EU - max. 24 pcs. for tray type FU - max. 48 pcs.
8854556	Clamps for Erlenmeyer flasks 250 ml required / possible max. number of clamps for tray type EU - max. 17 pcs. for tray type FU - max. 39 pcs.

Cat.-no.	Features and Specifications
8854572	Clamps for Erlenmeyer flasks 500 ml required / possible max. number of clamps for tray type EU - max. 12 pcs. for tray type FU - max. 26 pcs.
8854599	Clamps for Erlenmeyer flasks 1.000 ml required / possible max. number of clamps for tray type EU - max. 8 pcs. for tray type FU - max. 17 pcs.
8854610	Clamps for Erlenmeyer flasks 2.000 ml required / possible max. number of clamps for tray type EU - max. 4 pcs. for tray type FU - max. 9 pcs.
8854629	Clamps for Erlenmeyer flasks 3.000 ml, required / possible max. number of clamps for tray type EU - max. 4 pcs. for tray type FU - max. 8 pcs.
8854637	Clamps for Erlenmeyer flasks 5.000 ml required / possible max. number of clamps for tray type EU - max. 2 pcs. for tray type FU - max. 6 pcs.
	Clamps for Fernbach flasks (optional equipment of the universal trays)
8854564	Clamps for Fernbach flasks 450 ml required / possible max. number of clamps for tray type EU - max. 6 pcs. for tray type FU - max. 15 pcs.
8854600	Clamps for Fernbach flasks 1.800 ml required / possible max. number of clamps for tray type EU - max. 1 pc. for tray type FU - max. 6 pcs.
8854640	Clamps for Fernbach flasks 2.800 ml required / possible max. number of clamps for tray type EU - max. 1 pc. for tray type FU - max. 6 pcs.

Cat.-no.	Features and Specifications
----------	-----------------------------

**Plastic clamps for Erlenmeyer flasks
(optional equipment of the universal trays)**

8854700	Plastic clamps for Erlenmeyer flasks 100 ml required / possible max. number of clamps for tray type EU - max. 20 pcs. for tray type FU - max. 58 pcs.
---------	--

8854711	Plastic clamps for Erlenmeyer flasks 250 ml required / possible max. number of clamps for tray type EU - max. 20 pcs. for tray type FU - max. 40 pcs.
---------	--


8854722	Plastic clamps for Erlenmeyer flasks 500 ml required / possible max. number of clamps for tray type EU - max. 16 pcs. for tray type FU - max. 26 pcs.
---------	--

8854733	Plastic clamps for Erlenmeyer flasks 1.000 ml required / possible max. number of clamps for tray type EU - max. 9 pcs. for tray type FU - max. 15 pcs.
---------	---

Note: Plastic clamps are ideal for use with Erlenmeyer flasks with tubing connection. See below!

Shaking flask GL Erlenmeyer type,
with baffles for increased turbulence

Fig. 40:
Shaking flask Erlenmeyer type; with baffles in the concave side-wall section and connector GL14 for tubing connection

shaking flask made of DURAN[®]-glass, Erlenmeyer-type; straight neckpiece for metal caps; neckpiece of Ø = 38 mm; flask body with three baffles at Δ 120deg with GL 14 - connector incl. plastic cap and hose connectors for tubings of size 4 x 7 mm

8861064	size 300 ml, max. Ø = 87 mm, H = 161 mm; set of 10 pcs.
---------	---

8861072	size 500 ml, max. Ø = 105 mm, H = 183 mm; set of 10 pcs.
---------	--

8861080	size 1.000 ml, max. Ø = 131 mm, H = 232 mm; set of 10 pcs
---------	---


Accessories for baffles shaking flasks

8861099	cap made of aluminium, set of 10 pcs.
---------	---------------------------------------

8861102	stainless steel cap, set of 10 pcs.
---------	-------------------------------------

**Shaking flask Erlenmeyer type,
with baffles for increased turbulence**

Fig. 41:
Shaking flask Erlenmeyer
type;
with baffles in the
concave side-wall section


shaking flask made of DURAN®-glass, Erlenmeyer-type;
straight neckpiece for metal caps; neckpiece of Ø = 38 mm or version with
cellucotton stoppers
flask body with three baffles at Δ 120deg

- 8861005 size 300 ml, max. Ø = 87 mm, H = 161 mm; set of 10 pcs.
 - 8861013 size 500 ml, max. Ø = 105 mm, H = 183 mm; set of 10 pcs.
 - 8861021 size 1.000 ml, max. Ø = 131 mm, H = 232 mm; set of 10 pcs.
 - 8861022 size 2.000 ml, max Ø = 166 mm, H = 305 mm; set of 10 pcs.
 - 8860998 size 500 ml, version for cellucotton stoppers, etc.
max Ø = 105 mm, H = 180 mm; set of 10 pcs.
- Other versions on request

Universal Mounting Set
Basic element type B,
each consisting of 2 side frames and 4 rods


Fig. 42:
Basic element type B of
universal mounting set

- 8854238 type B-2 for trays EU
- 8854246 type B-3 for tray FU

Universal clamping rod type U


Fig. 43:
Universal clamping rod

8854254


type U for basic elements B-2 and B-3

Fig. 44:
Mounting system with
universal clamping rods

Mounting set for separating funnels type S

Fig. 45:
Mounting set for
separating funnels type S

8854262

type S-1 for 5 pcs. separating funnels 50 or 100 ml; set includes clamps, retaining springs and a universal clamping rod U-2

8854270

type S-2 for 3 pcs. separating funnels 250, 500 or 1.000 ml; set includes clamps, retaining springs and a universal clamping rod U-2

Max. number of mounting sets for tray type is

1 mounting set S 1 or S 2 for tray type EU with set B-2

2 sets of mountings S 1 or S 2 for tray type FU with set B-3

12.3 Safety Instructions Laboratory Shakers and Incubation Shaking Cabinets

Safety Instructions Laboratory Shakers and Incubation Shaking Cabinets see appendix)

12.4 EG Declaration of Conformity

(EG Declaration of Conformity see appendix)

12.5 Declaration of Decontamination

(Declaration of Decontamination see appendix)

12.6 Information and Instructions on Disposal and Repairs

(Information and Instructions on Disposal and Repairs see appendix)

12.7 Dimension sheet CERTOMAT® Tplus

(Dimension sheet CERTOMAT® Tplus see appendix)

12.8 Prospectus CERTOMAT® Tplus

(Prospectus CERTOMAT® Tplus see appendix)

Sartorius Stedim Biotech GmbH
August-Spindler-Straße 11
37079 Göttingen, Germany

Telephon +49.551.308.0
Fax +49.551.308.3289
www.sartorius-stedim.com

Copyright by
Sartorius Stedim Biotech GmbH,
Göttingen, Germany.
All rights reserved. No part
of this publication may be
reprinted or translated in any
form or by any means without
the prior written permission of
Sartorius Stedim Biotech GmbH.

The status of the information,
specifications and illustrations in
this manual is indicated by the
date given below.
Sartorius Stedim Biotech GmbH
reserves the right to make
changes to the technology,
features, specifications, and
design of the equipment
without notice.

Status:
July 2008
Sartorius Stedim Biotech GmbH,
Göttingen, Germany

Printed in Germany on paper that has
been bleached without any use of chlorine
W4A000 · KT
Publication No.: SE-6050-t08091
Order No.: 85032-536-02

Safety Instructions

Laboratory Shakers and Incubation Shaking Cabinets

Content

1	Introduction	3
2	Transport and Installation	3
2.1	Transporting Equipment	3
2.2	Requirements at the site of installation	3
3	Operating Instructions	4
3.1	Working with Media	4
3.2	Initial Startup and Normal Operation	4
4	Cleaning, Maintenance, and Service	5
4.1	Cleaning	5
4.2	Maintenance and Service	5
4.3	Returning Defective Units	5
4.4	Relocating equipment	5

1 Introduction

These Safety Instructions contain our recommendations for precautions in handling CERTOMAT® laboratory shakers and incubation shaking cabinets from the product program of the company Sartorius Stedim Biotech GmbH, especially for

- transport and installation
- equipping and initial startup
- operation
- cleaning, maintenance, and service

These safety instructions are non-transferable and meant exclusively for the equipment delivered. They are not valid for other equipment. Additional, specific safety regulations (i.e., legally or otherwise mandatory rules) are not dealt with by these instructions. Please consult your local authorities for further information.

Possible risks and dangers are marked with this symbol and highlighted like this paragraph. Ignoring these warnings may cause damage to the instrument or lead to other material or personal damage.

Steps to be carried out with special care or special aspects or issues to be considered are marked with safety instructions formatted as this paragraph.

[→ ..] Identifies references to contents of these instructions or other documents. Labels of illustrations, sections, or documents are in parentheses.

2 Transport and Installation

2.1 Transporting Equipment

The shakers and especially the incubation shaking cabinets are heavy. Use only suitable transportation and lifting tools when moving or setting up units.

Make sure clearances are sufficient along the transport route and that the combined weight of lifting and transportation tools does not exceed the maximum load permitted for the floor or ground.

Use only qualified personnel to transport the equipment. Make sure personnel is protected from injury during transport.

Do not remove any transport protection before the equipment has been transported to the site of installation.

2.2 Requirements at the site of installation

1. The floor or the laboratory bench where the equipment is to be installed must be able to support the weight of the fully equipped units and systems.
2. The installation place should be sufficiently dimensioned and feature a non-slip surface.
3. The equipment has to be easily accessible for operation, maintenance, and service tasks.
4. The shakers and the incubation shaking cabinets may transmit vibrations to the installation surface especially when carrying a heavy load. For this reason, install the shakers so that surrounding equipment is not affected.
5. The power supply must match the rating indicated on the type plate and must be grounded correctly.

Use a level to align the equipment after setting it up!

Before stacking equipment (max. 3 units), make sure that the bottom most device is properly aligned horizontally!

3 Operating Instructions

1. Make sure, that only authorized people have access to the working area.
2. Personnel has to be instructed about applicable safety instructions and other legally or otherwise compulsory regulations.

3.1 Working with Media

The utilized media and substances can be hazardous and may pose dangers that are specific to the substances and the processes and thus cannot be described here in detail. Your company should issue corresponding safety regulations and carefully instruct the operating personnel.

Recommended basic protection equipment includes such items as suitable working clothes, gloves, goggles, and respiratory equipment, if necessary.

3.2 Initial Startup and Normal Operation

1. Use only units, equipment, accessories, and spare parts released by the company Sartorius Stedim Biotech GmbH for use with the shaker.
2. Check all components for damage, especially glass components such as jars. Any damaged components must not be used.
3. Never use the shaker without the tray. Switch the unit off before installing/uninstalling trays, mounting sets, or additional accessories.
4. When mounting the trays, leave enough space between the unit and other objects, especially in case of trays extending over or beyond the shaker table or bench. Make sure that trays are properly installed and fastened.
5. Fasten the tray mountings and accessories properly. Distribute loads evenly on the tray when loading the shaker with the shaking vessels (Erlenmeyer flasks, bottles, etc.).
6. In case of incubation shaking cabinets, make sure the tray is firmly positioned; close the doors before starting the shaking process.
7. Before starting the shaker and during normal operation, check that the vessels are firmly positioned in their holders. The vessels must not become loose, rattle, or touch one another or the walls of the shaker. Switch the unit off immediately if vessels become loose in their holders.
8. Operate the equipment only within the permitted range of operating parameters. Comply with the equipment instructions and the technical data.

Head the lateral movement of the trays and the resulting impact or crushing hazards, especially with large trays, in case of low clearances to other equipment or objects or large amplitudes and heavy loads.

Never put fingers between tray support and frame when the unit is active. The space between these components is narrow due to constructional constraints. Risk of injury!

9. After switching off the shaker, wait until the shaker support has come to a complete standstill before removing or adding vessels.

4 Cleaning, Maintenance, and Service

Regular cleaning and maintenance ensures the equipment's functionality and safety.

-  Comply with the legal safety regulations applicable to the field of application concerning the handling of microbial material and contaminated equipment.
-  Prior to any cleaning or maintenance work, disinfect or sterilize any biologically contaminated components according to the safety regulations.
-  Prior to any cleaning or maintenance work, switch the unit off and pull the power cable to disconnect the unit from the mains power supply. Make sure that the unit cannot be restarted or activated accidentally.

Servicing and repairing drive systems, replacing drive belts, or working on electrical/electronic equipment should only be performed by qualified service personnel.

4.1 Cleaning

1. Clean the surfaces of the units with standard laboratory cleaners. Do not use corrosive agents (e.g., chlorous agents).
2. Broken glass or spilled liquids should be removed as quickly as possible.
3. For incubator shakers special hints are given in the operating manual.

4.2 Maintenance and Service

Maintenance and troubleshooting by the user should be restricted to the following tasks:

- Checking the voltage supply in the lab and the equipment connections.
- Checking and replacing defective fuses.
- For incubation shaking cabinets: Replacing defective light fittings.

If problems or malfunctions cannot be solved or remedied, please contact your Sartorius AG representative or contact us directly:

Sartorius AG
Servicezentrum Nord
Weender Landstraße 94-108
37075 Göttingen, Germany
Telephone +49.551.308.3729 / 3740 / 3741
Fax +49.551.308.3730

A written cleaning and decontamination certificate has to be provided when requesting service technicians of the company Sartorius AG to perform maintenance tasks.

4.3 Returning Defective Units

Defective units or components can be returned to the responsible Sartorius AG representative or directly to Sartorius AG.

1. The equipment has to be clean, hygienically flawless, and carefully packaged.
2. Contaminated parts have to be cleaned, disinfected, or sterilized according to the respective and applicable safety guidelines. Compliance with the guidelines has to be documented with a cleaning and decontamination certificate to be enclosed with the shipment.
3. Negotiate the return shipment with the Service Dept. before shipping any equipment.

4.4 Relocating equipment

1. Comply with the described cleaning and maintenance measures described to avoid biological or chemical hazards.
2. Disconnect the equipment from the mains power supply. Remove vessels and additional equipment or components inside of the unit or connected to the unit.
3. If applicable, reinsert the transport protection. Read and comply with the additional notes under "Transporting Equipment" (2.1).

Sartorius Stedim Biotech GmbH
August-Spindler-Straße 11
37079 Göttingen, Germany

Telephone +49.551.308.0
Fax +49.551.308.3289
www.sartorius-stedim.com

Copyright by
Sartorius Stedim Biotech GmbH,
Göttingen, Germany.
All rights reserved. No part
of this publication may be
reprinted or translated in any
form or by any means without
the prior written permission of
Sartorius Stedim Biotech GmbH.

The status of the information,
specifications and illustrations in
this manual is indicated by the
date given below.
Sartorius Stedim Biotech GmbH
reserves the right to make
changes to the technology,
features, specifications, and
design of the equipment
without notice.

Status:
February 2008,
Sartorius Stedim Biotech GmbH,
Göttingen, Germany

Printed in Germany on paper that has
been bleached without the use of chlorine
W4A000 · KT
Publication No.: SB6043az08022
Order No.: 85030-520-67

CE EC Declaration of Conformity

according to EU Machinery Directive 98/37/EC,
Appendix II A

Company	Sartorius Stedim Biotech GmbH	
Address	August-Spindler-Strasse 11 37079 Goettingen; Germany Telefon +49.551.308.0, Fax +49.551.308.3289 www.sartorius-stedim.com	
	<p>We herewith declare that the device described below fulfils the relevant fundamental safety requirements and health regulations specified by the appropriate EU-Directive, with respect to its design and construction and to the version as commercialized.</p> <p>This declaration becomes legally invalid if modifications are performed on the device which have not been certified by Sartorius Stedim Biotech GmbH.</p>	
Designation of the device	CERTOMAT® T plus	
Model, version	Incubation Shaking Cabinet	
Cat.-No.	8865809, 8865825, 8865906, 8865922	
Relevant directives of the EC	Directive 98/37/EC	Machinery Directive
	Directive 2004/108/EC	Electromagnetic Compatibility
	Directive 2006/95/EC	Electrical equipment designed for use within certain voltage limits
Applied standards	EN ISO 12100-1:2003, EN ISO 12100-2:2003 EN 61326-1:2006 EN 61010-2-010:2001	
Date and Signature	20.08.08 <i>I.V. P. B.</i> <i>S. Gerighausen</i>	
Function of Signatory	Dr. Reinhard Baumfalk Vice President Operations Instruments	Dr. Susanne Gerighausen Quality Management Representative

Declaration about decontamination and cleaning of equipment and components

When returning equipment or components, please describe on page 2 of this form the problem(s) or fault(s) you have found. Please also indicate the remedial actions you require.

To protect our personnel, we require all equipment or components be free of biological, chemical, or radioisotopic contaminants. We will only accept such equipment or components when:

- the equipment or components have been adequately cleaned and decontaminated.
- this declaring document has been completed, signed and returned by an authorized person.

Please help us in assuring a safe, hazard-free work environment.

A. Description of the Equipment or Component(s)

Description / Cat. No.	
Serial no.	
No. of invoice/delivery note	
Date of delivery	

B. Contamination / Cleaning

Attention: Please specify exactly the biological, chemical, or radioisotopic contaminant	Attention: Please describe the cleaning and decontamination procedure/method
The equipment was contaminated with	and it has been cleaned and decontaminated by

F-13-01 Version 01	Page 2 / 4 Effective Date 01.01.2008	<h2>Return of material</h2>	
------------------------------	---	-----------------------------	--

C. Legally binding declaration

I (we) certify that all information given in this form is correct and complete.
 The equipment and components have been adequately decontaminated and cleaned according to the legal requirements.
 No chemical or biological or radioisotopic risks remain that can endanger exposed persons' safety or health.

Company / Institute	
Address / Country	
Tel. / Fax (with area code)	
Name of the authorized person	
Position	
Signature / Date	

D. Reason for return

wrong delivery
 exchange
 repair
 modification
 disposal
 other

E. Please describe the problem(s) or fault(s) you have found (for repair) and/or indicate the remedial actions you require

F-13-01	Page 3 / 4	Return of material	
Version 01	Effective Date 01.01.2008		

F. Reserved for Sartorius-Service Center

Notes

Please pack the unit properly and send it freight paid to your local service supplier or directly to Sartorius AG Servicezentrum Nord, Germany.

Sartorius AG
Servicezentrum Nord
Weender Landstraße 94-108
37075 Göttingen, Germany
Telephone +49.551.308.3729 / 3740 / 3741
Fax +49.551.308.3730

F-13-01 Version 01	Page 4 / 4 Effective Date 01.01.2008	Return of material	 sartorius stedim biotech
--------------------------	---	---------------------------	---

Sartorius Stedim Biotech GmbH
August-Spindler-Straße 11
37079 Göttingen, Germany

Telephone +49.551.308.0
Telefax +49.551.308.3289
www.sartorius-stedim.com

Copyright by
Sartorius Stedim Biotech GmbH,
Göttingen, Germany.
All rights reserved. No part
of this publication may be
reprinted or translated in any
form or by any means without
the prior written permission
of Sartorius Stedim Biotech
GmbH.

The status of the information,
specifications and illustrations in
this manual is indicated by the
date given below. Sartorius Stedim
Biotech GmbH reserves the right
to make changes to the technology,
features, specifications, and
design of the equipment
without notice.


Status:
February 2008,
Sartorius Stedim Biotech GmbH,
Göttingen, Germany

Information and Instructions on Disposal and Repairs

Packaging that is no longer required must be disposed of at the local waste disposal facility. The packaging is made of environmentally friendly materials that can be used as secondary raw materials.

The equipment, including accessories and batteries, does not belong in your regular household waste. The EU legislation requires its Member States to collect electrical and electronic equipment and disposed of it separately from other unsorted municipal waste with the aim of recycling it.

In Germany and many other countries, Sartorius AG takes care of the return and legally compliant disposal of its electrical and electronic equipment on its own. These products may not be placed with the household waste or brought to collection centers run by local public disposal operations – not even by small commercial operators.

For disposal in Germany and in the other Member States of the European Economic Area (EEA), please contact our service technicians on location or our Service Center in Goettingen, Germany:


Sartorius AG
Service Center
Weender Landstrasse 94-108
37075 Goettingen, Germany

In countries that are not members of the European Economic Area (EEA) or where no Sartorius affiliates, subsidiaries, dealers or distributors are located, please contact your local authorities or a commercial disposal operator.


Prior to disposal and/or scrapping of the equipment, any batteries should be removed and disposed of in local collection boxes.

Sartorius AG, its affiliates, subsidiaries, dealers and distributors will not take back equipment contaminated with hazardous materials (ABC contamination) – either for repair or disposal. Please refer to the accompanying leaflet/manual or visit our Internet website (www.sartorius.com) for comprehensive information that includes our service addresses to contact if you plan to send your equipment in for repairs or proper disposal.

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CERTOMAT Tplus: ca. 200kg

Untergestell Tplus: ca. kg

The reproduction, distribution and utilization of this document as well as the
 creation of copies in any form or by any means, including electronic, is
 prohibited. All rights reserved in the event of a patent, utility model or design.
 Änderungen nur über CAD!

Index	Field	Changes/Änderung

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Sartorius Stedim Singapore Pte. Ltd.
Sartorius Stedim Australia Pty. Ltd.